

CATALOG

2023 EDITION



ANESTHESIA & AIRWAY
GASTROENTEROLOGY
LABORATORY PRODUCTS
ORTHOPEDIC PRODUCTS
RENAL CARE

RESPIRATORY CARE STERILIZATION SURGERY & OPERATION UROLOGICAL PRODUCTS









Who we are

Inmed Corporation is the country's leading and trusted master distributor of quality local and imported medical products and supplies. With the primary focus of providing solid marketing support, Inmed Corporation has established a dependable structure offering significant value to its resellers and dealers in the Philippines and neighboring Asia Pacific Countries.

Committed to help both local and multinational medical partner companies in the Philippines to achieve their target goals. Inmed Corporation continues to work hand-in-hand with them to extend and further expand market reach. The company utilizes its financial strength to purchase high volume supplies and offer its resellers savings at the best attainable cost.

Managed by a team of top professionals in the industry with combined experience of over 60 years in medical products, sales, marketing and distribution, Inmed Corporation boasts of seasoned round-table negotiations, financial packaging and the best business practices benchmarked on international standards providing economical and fundamentally sound solutions to your business needs.

Our Mission

Our mission is to provide quality local and imported medical products at competitive prices and offer professional customer-oriented service.

Our Vision

Our vision is to be the Philippines' leading master distributor of quality medical products — being the preferred supplier of all resellers in the industry and recognized as a strong partner in sustainable healthcare in the nation.

ISO Certified

INMED Corporation is ISO 9000:2015 certified by TUV, compliant with the International Organization for Standardization based in Geneva, Switzerland. It is our commitment to provide our customers with only high quality products, unparalleled service excellence, and reliable timely delivery of all customer orders. Customer satisfaction is our guarantee.





TABLE OF CONTENTS

2 & 3 Port Manifold	62	Ethylene Oxide (EO) Indicator Strip	46
3 Port Manifold	62	Feeding Tube	14
Abdominal Poole Drain Suction Set	57	Feeding Tube with X-Ray Line	14
Absorbable Gelatin Sponge	54	Foley Balloon Catheter (Silicone Coated)	67
Absorbent Surgical Gauze Roll	53	Frazier Tip with Suction Connecting Tube	57
Aerosol Mask	34	Gastric Sump Tube	17
Amber Rubber Tubing	22	Gastrostomy Feeding Tube	19
Anesthesia Breathing Circuit	11	Gauze Ball with X-Ray Line (Cherry Sponge)	53
Anesthesia Mask	10	Gauze Ball with X-Ray Line (Peanut Sponge)	53
Anesthesia Mask Head Strap	10	Gauze Swab/Sponge (Sterile & Non-Sterile)	53
Angiography Drape	51	Gigli Saw Wire	62
Angiography Pack	52	Guedel Airway	6
Asepto Bulb Syringe	18	Green Top Heparin - 2ml	26
Autoclave Tape	44	Heat and Moisture Exchanger (HME) Filter	40
AV Fistula with Back-Eye Needle	31	Hemocytometer Cover Glass	23
Bacterial Viral Filter (BVF)	40	Hemodialysis Catheter Kit	30
Blood Bag	31	High Flow Nasal Cannula	33
Blood Lancet	21	High-Concentration Oxygen Mask with Tubing	34
Blood Transfusion Set	31	High-Flow Humidification System	42
Blood Tubing Set	30	Hollow Fiber Dialyzer High Flux	30
Blue Top Sodium Citrate - 2ml	27	Hollow Fiber Dialyzer Low Flux	30
Bone Wax	54	Hydrophobic Suction Filter	62
Bouffant Doctor's Cap	50	Incentive Respiratory Exerciser	41
Bouffant Nurse Cap	50	Inoculating Loops	22
Bowie and Dick Test Pack	44	Intubating Stylet	9
Cathy Closed Suction System	38	Laparotomy Pack	51
		Laparotomy Sponge	54
Cautery Tip Cleaner	55 52	Laryngeal Mask Airway	9
Cesarean (CS) Pack	60	Latex Foley Balloon Catheter	66
Chest Drainage		Lavender Top EDTA - 3ml and 5ml	26
Chest Drainage Thoracic Catheter	60	M Lab Sealing Film	23
Closed Suction Catheter System	37	Maintenance Care Pack	39
Closed System Wound Suction Kit	59	Male External Condom Catheter	67
Corrugated Aerosol Tubing	40	Manual PVC Resuscitator	36
Cover Glass	23	MDI Adapter	40
Crepe Paper Wrapping Sheet	46	Medium Concentration Oxygen Mask with Tubing	34
Disposable Apron	21	Micro Slides	23
Disposable Lab Gown	21	Mounting Medium	24
Disposable Skin Staple Remover	56	Mouthpiece with Elastic Band	17
Disposable Surgical Gown	49	Mucus Extractor	38
Double Lumen Endobronchial Tube	9	Multi-Draw Needle - G-21 x 1 1/2"	27
Duodenal Tube	15	MIC-KEY Low-Profile Gastrostomy Feeding Tube	19
Duodenal Tube with X-Ray Line	16	Multi-Vent Mask	34
Electrosurgical Cautery Pencil	55	Nasal Oxygen Cannula	33
Electrosurgical Pad	55	Nasopharyngeal Airway	12
Endotracheal Tube - Cuffed & Uncuffed	7	National Scrub Brush	63
Endotracheal Tube Holder	6	Nebulizer Kit with Mask	35
Endotracheal Tube Introducer	10	Nebulizer Kit with Mouthpiece	35
Endotracheal Tube with Subglottic Port	8	Nelaton Catheter	65
Enteral Feeding Container	16	Obstetrics (OB) Pack	51
Enteral Replacement Tubing Set	16	Ophthalmic Pack	52

TABLE OF CONTENTS

Orthopedic Cast Padding	28	Toomey Syringe	18
Oxygen Connecting Tube	33	Trachcare Closed Suction Catheter	38
Peak Flow Meter	41	Tracheostomy Mask	39
Pediatric Urine Collector	67	Transfer Pipette	22
PEG Feeding Tube Kit	14	Tri-Spring Wound Suction Kit with Trocar	58
Penrose Drain Tubing	59	T-Tube Drain	59
Petri Dish	25	Tube Needle Holder	21
Pipette Tips	22	Turbo Cleaning Closed Suction System	37
Plasma Indicator Strip	46	Tyvek Sterilization Flat Reel	45
Plaster of Paris	28	Umbilical Cord Clamp	56
Pre-filled Humidifier 1000ml	35	Umbilical Cord Clamp Clipper	56
Pre-Filled Humidifier with Adaptor	35	Unigard Urine Drainage Bag with non-return valve	68
Preformed Endotracheal Tube - Oral & Nasal	8	Universal Funnel Tubing	58
Pre-Loaded Endotracheal Tube	8	Universal Pack	52
Pre-Operational Disposable Razor	55	Unometer Safeti Plus	69
Rebreather Bag	11	Urine Bedside Bag	68
Rectal Catheter	17	Urine Leg Bag	68
Red Nelaton Catheter	66	Urine Meter	68
Red Top Plain - 3ml and 5ml	26	Urine Test Strips	25
Reinforced Endotracheal Tube	9	Ventilator Circuit	41
Safeliner	61	Vessel Loop	61
Salem Sump	18	Video Laryngoscope	10
SARS-CoV-2 Antigen Test Kit	25	Yankauer Suction Set	57
Self-Sealing Sterilization Pouch	44	Yankauer Suction Set	57
Shoe Cover (Non-Skid)	50 15	Yellow Top Clot Activator with Gel - 5ml	26
Silicone Duodenal Tube	15		
Silicone Foley Balloon Catheter Skin Stapler - W35	66 56		
·	36 47		
SMS Wrapping Sheet Soda Lime	47 11		
Soft Touch Tracheostomy Tube Holder	39		
Solidifier Gel	61		
Specimen Container	24		
Spinal Needle Quincke Point	12		
Spirometer Filter with Nose Clip	39		
Steam Autoclave	44		
Steam Indicator Strip	47		
Sterile Surgical Blade	54		
Sterile Surgical Gloves	50		
Sterilization Flat Reel	45		
Sterilization Flut Reel	45		
Suction Catheter	36		
Suction Catheter with Finger Tip Control	37		
Suction Connecting Tube	58		
Suction Connecting Tube	58		
Surgical Gown	49		
Surgical Needles	63		
Thoracic Catheter	60		
Tissue Embedding Cassette	24		

24

Tissue Embedding Medium



Guedel Airway

Sedasenz Guedel Airway is a medical device used as an airway adjunct to maintain the patient's airway open.

FEATURES AND BENEFITS:

- ► Made from medical grade Polyethylene (PE) material
- ► Soft rounded edges to minimize trauma from the patient's tissue
- ► Anatomically correct shape
- Color-coded based on sizes

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available in Size: 00, 0, 1, 2, 3, 4 & 5	Ea	50	500	05709



simplex

Guedel Airway

Simplex Guedel Airway provides free flow of air while protecting the airway from occlusion. It secures an unobstructed oropharyngeal airway.

FEATURES AND BENEFITS:

- ► Anatomically correct shape follows the contour of the palate
- ► Soft rounded edge prevents tissue damage
- ► Color-coded for easy size identification
- ► Bite block provides access for suction catheters

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available in Size: 00, 0, 1, 2, 3 & 4	Ea	50	500	5779



Sedasenz[®]

Endotracheal Tube Holder

Sedasenz Endotracheal Tube Holder is a stabilization device that secures the patient's endotracheal tube in place, avoiding complicated kinks, tube dislodgement and time-consuming fixation with plaster tapes.

- ► Bite block protects patient's teeth and gums allowing continuous airflow.
- ► Push-screw mechanism allows easy tighten or release of the tube clamp.
- ► Fits any size of endotracheal tube.
- ► Soft neck-tie with Velcro strap that allows easy adjustment.
- ► With inner foam padding for patient's comfort

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-PM-ETHD01	Sedasenz Endotracheal Tube Holder	Ea	25	150	_



Endotracheal Tube Cuffed & Uncuffed

Sedasenz Endotracheal Tube is inserted either through the nose or through the mouth and into the trachea. It is connected to the ventilation device to maintain patient's airway.

FEATURES AND BENEFITS:

- With 15mm standard connector to ensure compatibility with circuit connectors
- Radiopaque line which is visible during X-ray for confirmation of placement
- With HVLP cuff to atraumatically and securely seal the tracheal airway, preventing leakage and gastric aspiration




DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Without Cuff 2.0 to 4.5	Ea	10	100	05525
With Cuff 5.0 to 9.0	Ea	10	100	05525

simplex

Endotracheal Tube Cuffed & Uncuffed

Simplex Endotracheal Tubes is designed to provide artificial airway to patients by allowing positive pressure management and tracheobronchial suctioning. Available in wide range of variants: Pediatric (with and without cuff) and Adult (with and without cuff).

- ► Thermosensitive material conforms to the patient's body temperature, reducing trauma to the patient's tissue
- ► Kink-resistant shaft allows continuous air flow
- ► Radiopaque line allows verification of tube placement through X-ray
- ► Responsive pilot balloon allows visual identification of the cuff's integrity
- ► Double line marking facilitates accurate cuff placement below the vocal cord
- High Volume Low Pressure (HVLP) Cuff that provides effective seal, reducing risk of vocal cord injury and tracheal damage
- ► Murphy eye at the tip for increase safety, ensuring adequate ventilation is maintained



DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Without Cuff 2.0 to 6.0	Ea	10	100	5777
With Cuff 3.5 to 8.5	Ea	10	100	5774

Preformed Endotracheal Tube Oral & Nasal

Sedasenz Preformed Endotracheal Tubes are molded in special shapes that permits good surgical access in the oro-nasal area.

- Preformed South/Oral provides good access for the ENT surgeons needing to work on the head, opthalmic or nasal area.
- Preformed North/Nasal provides good access to the mouth for dental work or for any mandibular operations.

FEATURES AND BENEFITS:

- ► Preformed curvature reduces the risk of kinking and tube obstruction
- Rounded tip with murphy eye is smoothly formed to avoid vocal cord damage and to ensure atraumatic and smooth intubation
- ► High-volume low-pressure (HVLP) cuff is designed to provide effective sealing & increase pressure distribution
- ► Double marking to guide in intubation to help ensure accurate placement
- ► Radiopaque line facilitates verification of tube positioning via X-ray

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Oral Curved Endo Tube w/ Cuff, 7.0 to 8.0	Ea	10	100	05525
Nasal Curved Endo Tube w/ Cuff, 7.0 to 8.0	Ea	10	100	05525



Nasal Preformed

Sedasenz[®]

Endotracheal Tube with Subglottic Port

Sedasenz Endotracheal Tube with Subglottic Port enables elimination of secretions above the cuff, designed to prevent ventilator associated pneumonia.

FEATURES AND BENEFITS:

- ► Subglottic port allows removal of subglottic secretions using a syringe or suction equipment
- ► Rounded bull-nose tip with murphy eye is smoothly formed to avoid vocal cord damage
- ► HVLP cuff provides effective sealing and increases pressure distribution
- ► Double markings to guide in intubation to help ensure accurate placement
- ► Radiopaque X-ray line facilities verification of tube placement

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available in 7.0 to 8.0	Ea	10	100	05525



Sedasenz[®]

Pre-Loaded Endotracheal Tube

Sedasenz Pre-loaded Endotracheal Tube is an endotracheal tube with stylet incorporated inside the tube, it is designed for fast set up to achieve successful emergency intubation.

- ► 15mm connector ensures compatibility with circuit connectors
- ► Radiopaque line for confirmation of the tube placement
- Stylet has low friction coating to facilitate smooth withdrawal
- HVLP cuff for prevention of air leakage & gastric aspiration
 With soft distal tip to reduce risk of mucosal damage

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Uncuffed, 2.5 to 4.5	Ea	10	100	05525
Cuffed, 5.0 to 8.5	Ea	10	100	05525



Reinforced Endotracheal Tube

FEATURES AND BENEFITS:

- ► Durable spiral wire embedded in the endotracheal tube provides effective resistance to kinking
- Equipped with murphy eye at the distal tip for continuous ventilation even when distal tip is occluded
- ► High-volume low pressure cuff, designed to provide effective sealing and increase pressure

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: 7.0 - 8.0	Ea	5	200	05525



Sedasenz[®]

Double Lumen Endobronchial Tube

FEATURES AND BENEFITS:

- ► Adapter enables a simple one-step connection to the breathing system
- Color-coded & printed pilot balloons with 1-way valves for clear indications of cuff inflation routes & status

- Clear distance markings assists accurate positioning of the tube
 High Volume Low Pressure Cuff minimize trauma to tracheal mucosa while achieving an excellent seal
 Manufactured from clear, thermosensitive PVC material which softens at patient's body temperature, preventing trauma to delicate tracheal mucosa
- Larger lumen allows use of larger fiberoptic bronchospes and suction catheters

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available in FR-35 and FR-37	Ea	1	20	05110



Sedasenz[®]

Laryngeal Mask Airway

FEATURES AND BENEFITS:

- ► Made from medical grade silicone, enhancing patient comfort and superior biocompatibility.
- ► Intuitive design allows easy insertion.
- Soft cuff provides an airtight seal, allowing a secure airway.
- ► Resposive pilot balloon for easy cuff inflation

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Size available: 1 to 5	Ea	5	20	_



Sedasenz[®]

Intubating Stylet

- ► Made from malleable aluminum, covered with plastic sheath with a lubricated silk-slip surface
- Helps reduce friction between stylet and tracheal tube for easy insertion and withdrawal
- Sheath tip is rounded to reduce the risk of trauma if the stylet extends beyond the tip of the endotracheal tube

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available size: Fr-6 to Fr14	Ea	20	400	5775




Endotracheal Tube Introducer

Sedasenz Endotracheal Tube Introducer / Bougie is an adjunct for difficu allows atraumatic tube replacement.

FEATURES AND BENEFITS:

- ► Atraumatic rounded tip to reduce potential patient trauma
- ► Low friction surface for easy insertion and operation
- Coude tip facilitates correct placement by providing accurate feel of the tracheal wall c-rings
- ► Flexible to allow easy removal once the ET tube is in place

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-PM-SETB15	Sedasenz Endotracheal Tube Introducer, Fr-15	Ea	10	200	06106



MC

Video Laryngoscope

MC Video Laryngoscope is a convenient and durable video display tandem with disposable blades, engineered to provide confidence, precise, quick and successful endotrac heal tube intubation.

FEATURES AND BENEFITS:

- ► Responsive and easy to use 3.5" intuitive touchscreen with adjustable LED source and display luminosity
- Portable and waterproof
- HD image viewing with automatic focusing and anti-fog technology guarantees optimal visual capability
 Real time recording of video sequence and still image capture
- ► With comprehensive range of Macintosh blades for different types of patients
- ► Suitable for infant, children, adult, obese as well as difficult-to-intubate patients

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
MC Video Laryngoscope Main Unit	Unit	1	1	_
Laryngoscope Blade Size #1 to 4	Ea	12	120	_



Sedasenz[®]

Anesthesia Mask

Sedasenz Anesthesia Mask features an anatomically designed mask that provides non-invasive administration of anesthetic gas to a patient through inhalation during the induction of general anesthesia.

FEATURES AND BENEFITS:

- ► Smooth material allows tight seal while providing patient comfort
- Clear mask with inflatable cushion preventing pressure points on the patient's face
- ► Luer slip inflation valve for cushion adjustment
- ► Single-use

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Sedasenz Anesthesia Mask #2 and #5	Ea	1	100	05704



Sedasenz[®]

Anesthesia Mask Head Strap

Sedasenz Anesthesia Mask Head Strap is a mask harness for the anesthesia mask, designed to secure the anesthesia mask onto the patient's head.

- Four straps with holes that attach securely to the hook rings of the anesthesia mask.
- ► Fully adjustable mask harness
- ► Ideal for adult size of anesthesia mask

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-AMSHS1	Sedasenz Anesthesia Mask Head Strap, Adult	Fa	1	_	_



Anesthesia Breathing Circuit

Sedasenz Anesthesia Breathing Circuit is a complete assemble of components which connects the patient's airway to the anesthetic machine, creating an artificial atmosphere from and into which the patient breaths.

FEATURES AND BENEFITS:

- ► Universal connectors and flexible circuits for end user's efficiency
- ► Universal elbow connector with gas sampling port
- ► With extra limb for added connection
- ► 1 piece reservoir bag to meet the peak inspiratory flow requirement



DESCRIPTION / SIZE	QTY/PACK	QTY/CTN	MDR
2-Limbs, Adult and Pedia	1	20	05575
3-Limbs, Adult and Pedia	1	20	05575

Sedasenz[®]

Rebreather Bag

Sedasenz Rebreather Bag is a collapsible gas reservoir which permits manual ventilation and it also acts as a visual indicator of the patient's respiratory rate.

FEATURES AND BENEFITS:

- ullet Hourglass-shaped contoured breathing bag efficiently delivers the required gas volume
- ► Thin and durable allows tactile sensitivity
- ► Non-slip grip, textured exterior surface
- ► Latex-free

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in: 1L, 2L and 3 L	Ea	1		05575



Sedasenz[®]

Soda Lime

Sedasenz Soda Lime is used to absorb carbon dioxide in a closed breathing environment to prevent carbon dioxide retention which could lead to carbon dioxide poisoning.

- A granular mixture of calcium hydroxide with sodium hydroxide and color indicator
- ► White color changes to purple
- ► Granule hardness is 99%
- ► 3% moisture absorption

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-ABSL45	Sedasenz Soda Lime, 10lbs. (4.5kg)	Container	1	2	MDRE
0-PM-ABSL38	Sedasenz Soda Lime, 38lbs. (17.24kg)	Pail	1	1	MDRE





Spinal Needle Quincke Point

Top Spinal Needle Quincke Point enables smooth penetration and efficient administration of anesthetic medication, designed with perfect angular fit to prevent coring and damage to the body tissues.

FEATURES AND BENEFITS:

- ► Key/Slot arrangement style of the cannula hub facilitates proper needle bevel orientation
- ► Translucent window hub for easy visualization of cerebrospinal fluid
- ► Color-coded spinal needle gauge for size identification
- ► Needle length is 3-1/2 inch

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: G18 to G27	Ea	20	200	2995





Respisenz®

Nasopharyngeal Airway

Respisenz Nasopharyngeal Airway is a soft, flexible, uncuffed tube placed through the nasal passages so that the distal tip rests in the nasopharynx. It is used to maintain the free passage of air to and from the lungs in patients with facial trauma or lockjaw or in nearly comatose patients. Ideal to use when patient's oral passage is not accessible, it also allows nasotracheal suctioning.

- ► Soft PVC material, allows gentle insertion
- ► Flexible for maximum patient comfort
- ► Trumpet design for secure placement, preventing slipping of tube through the nasal passage

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: 6.0mm, Fr-24 to Fr-32	Ea	10	100	06089





ΔVΔNOS

PEG Feeding Tube Kit

Avanos Percutaneous Endoscopic Gastrostomy (PEG) Feeding Tube Kit is designed for placement convenience. Comes in complete placement kit for pull placement method.

FEATURES AND BENEFITS:

- ► Medical grade silicone construction
- ► Multiple ports accomodate nutrition and medication delivery
- ► Collapsible internal retention bumper
- ► SECUR-LOK external retention ring

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-20 and Fr-24	Ea	1	1	7536



simplex

Feeding Tube

Simplex Feeding Tube is used through nasal introduction for the purpose of administration of enteral nutrition and medication or for aspiration of intestinal fluids for newborn and infants.

FEATURES AND BENEFITS:

- ► Frosted surface eliminates need for lubricant
- ► With soft rounded tip for atraumatic introduction
- ► Also available: Simplex Feeding Tube with X-ray line and calibration markings

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-5, Fr-8 and Fr-10	Ea	100	1000	01629





Feeding Tube with X-Ray Line

Unomedical Feeding Tube with X-Ray Line is used for enteral feeding to provide nutrition to patients who cannot obtain nutrition by mouth, unable to swallow safely or need nutritional supplementation. Used also for administration of medicine to premature babies or small infants.

- ► Made from 100% DEHP-Free material
- Frosted surface provides smooth insertion of the feeding tube, less traumatic to the patient's tissue
- ► Kink-resistant tubing
- ► Detectable x-ray line and tip ensures position control
- ► Uno Metric tube marks at each centimeter for clear identification of insertion depth
- ► Displaced eyes minimizes risk of dumping syndrome
- Cap on funnel for protection against contaminants

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: Infant Fr-5, Fr-8 and Fr010	Ea	100	1000	_
Available Size: Adult Fr-5, Fr-8 and Fr010	Ea	100	600	



Unimex®

Duodenal Tube

Unimex Duodenal Tube is used through nasal introduction for direct to stomach administration of enteral nutrition or for aspiration of intestinal fluids.

FEATURES AND BENEFITS:

- ► Made from medical grade PVC
- ► With side eyes to reduce risk of blockage and for maximum dissemination
- ► Color-coded funnel for easy identification of size
- ► Markings along the tube to determine the depth of insertion

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Availabe Size: Fr-10 to Fr-18	Fa	50	500	05527



Unimex⁶

Silicone Duodenal Tube

Unimex Silicone Duodenal Tube is used through nasogastric introduction for administration of enteral nutrition or for aspiration of intestinal fluids.

FEATURES AND BENEFITS:

- ► 100% silicone material
- ► Cap on funnel, prevents bacteria from entering the tube if not in use
- ► Markings on the tube to determine the depth of insertion
- ► With radiopaque line and tip
- ► Length: 125cm

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Availabe Size: Fr-14 to Fr-18	Ea	50	500	05527



simplex

Duodenal Tube

Simplex Duodenal Tube is inserted into the nose and ends into the stomach, used to provide nutritional support to patients who are unable to ingest and swallow food.

- ► Made from transparent kink-resistant material allowing easy observation
- ► Frosted PVC material, prevents tissue damage
- ► Large lateral eyes for effective enteral feeding
- ► Color-coded based on french size
- ► Also available; Simplex Duodenal Tube with Funnel Cover

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Availabe Size: Fr-10 to Fr-18	Ea	50	500	01624
Availabe Size: Fr-10 to Fr-18 (with Cover)	Fa	50	500	01624







Duodenal Tube with X-Ray Line

Unomedical Duodenal Tube with X-Ray Line is inserted through nasogastric introduction, used for administration of enteral nutrition or for aspiration of intestinal fluid.

FEATURES AND BENEFITS:

- ► Transparent material, facilitates observation of aspirates or food passage
- ► Markings along the tube to determine the depth of insertion
- ► Detectable x-ray line and tip ensures placement verification via x-ray
- ► 4 lateral eyes with open tip, facilitates effective drainage

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-12 to Fr-18	Ea	100	400	04259



Unimex[®]

Enteral Feeding Container

Unimex Enteral Feeding Container is used during gravity feeding for the delivery of nutritionally complete food for malnourished or at risk of malnutrition patients, critically ill patients or post operative patients with limited oral intake.

FEATURES AND BENEFITS:

- ► Rigid container, stable structure and non-collapsible
- Large top opening, minimizes formula spills and waste
- Calibrated with easy to read graduations up to 600cc with 10cc increment, for improved feeding accuracy
- ► Clear drip chamber for visual identification of content
- ► 95cm tubing for greater patient comfort
- ► With roller flow control to regulate feeding rate
- ► Non-toxic, non-pyrogenic, non-latex, no medical interaction and safe to use
- ► Built-in strong, dependable and universal hanger for gravity feeding

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-PM-UEFC60	Unimex Enteral Feeding Container, 600cc	Ea	1	50	MDRE



Unimex[®]

Enteral Replacement Tubing Set

Unimex Enteral Replacement Tubing Set is a replacement tubing for Unimex Enteral Feeding Container.

- ► Made from medical grade PVC
- ► Non-toxic, non-pyrogenic
- ► Clear tubing with regulator, 95cm in length

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-PM-UEFT60	Unimex Enteral Replacement Tubing Set	Ea	1	300	



medovations.

Gastric Sump Tube

Medovations Gastric Sump Tube provides safer, continuous and intermittent gastric suctioning.

FEATURES AND BENEFITS:

- ► Double Lumen Tube that allows irrigation and suctioning
- ► Funnel end comes with universal step connector
- ► With radiopaque stripe for placement verification



DESCRIPTION / SIZE UOM		QTY/BOX QTY/CTN		MDR
Availabe Size: Fr-12 to Fr-18	Ea	1	50	05232



Rectal Catheter

Unomedical Rectal Catheter is inserted via the rectum to drain waste, or to perform lavage or for the purpose of medication administration in the lower intestine.

FEATURES AND BENEFITS:

- ► Made from medical grade Polyvinyl Chloride, biocompatible to human tissue
- ► Frosted surface reduces surface friction
- ► Closed rounded distal end with two off-set lateral eyes permits even distribution of liquid
- ► Color-coded funnel which allows easy identification of size
- ► DEHP-free

DESCRIPTION / SIZE		QTY/BOX		MDR
Available Size: Fr-20 to Fr-28	Ea	100	600	04373





Mouthpiece with Elastic Band

Greatcare Mouthpiece with Elastic Band is used to provide a passage for endoscopic insertion while preventing costly bite damage to the scopes during endoscopic procedures.

- ► Novel design
- Groves in the middle prevents the endopiece from moving
- ► With tongue depressor feature
- ► Rounded edges to provide patient comfort
- ► Hygienic and convenient

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR	
1-UN-MPP000	Greatcare Mouthpiece with Elastic Band	Ea	50	400	MDRE	





Asepto Bulb Syringe

Simplex Asepto Bulb Syringe is designed to fit directly into a large lumen tubing, used for intraoperative irrigation, and it can also be used to deliver nutritional support to a patient with a feeding tube.

FEATURES AND BENEFITS:

- ► Large and elongated tip to facilitate feeding and water flushes
- ► Smooth tip to minimize injury during wound irrigation
- ► Soft and pliable bulb for easy control
- ► 60cc capacity

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
2-PP-SBAS00	Simplex Asepto Bulb Syringem 60cc	Ea	1	50	882505



simplex

Toomey Syringe

Simplex Toomey Syringe is conveniently designed for irrigating and evacuating of fluids during medical procedures, can also be used during enteral feeding.

FEATURES AND BENEFITS:

- ► With plain style-tip, long and tapered to fit inside the lumen of a catheter tube
- ► Has a thumb ring plunger for one-hand use
- ► 60cc capacity graduated barrel

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-SP-TS0060	Simplex Toomey Syringe, 60cc	Ea	1	50	753144



_simplex_Salem Sump

Simplex Salem Sump is a double lumen tube with large lumen for aspiration and a small lumen for venting and irrigation, used to prevent the accumulation of gastric secretions or air from the stomach.

- ► Distal end is coned with ten lateral eyes
- ► Color-coded venting and irrigation port for easy identification
- ► Radiopaque tube is marked at 45, 55, 65 and 75cm from the tip for accurate placement
- Proximal end is provided with universal funnel connector along with double taper connector for easy extension
- ► Smooth and soft distal tip and tubing prevents mucosal damage

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Availabe Size: Fr10 to Fr-18	Fa	50	200	_





Gastrostomy Feeding Tube

Simplex Gastrostomy Feeding Tube is intended to provide means of feeding/medication administration directly to the patient's stomach.

FEATURES AND BENEFITS:

- ► Made from 100% medical grade silicone construction
- ► Graduated shaft allows easy stoma depth measurement
- ► With external bolster which allows air circulation, while preventing pressure from the stoma
- ► Flexible balloon construction provides better anchorage to prevent tube dislodge.

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-14 to Fr-24	Box	1	10	06996



ΔVANOS

Gastrostomy Feeding Tube

Avanos Gastrostomy Feeding Tube is indicated for patients who require enteral feeding but are unable to tolerate oral feeding. Gastrostomy tubes are inserted directly into the stomach and well suited for long term enteral feeding.

FEATURES AND BENEFITS:

High-clarity medical grade silicone construction for visibility of flow Ventilated external retention ring allows air circulation around the stoma site Clearly marked feeding, medication and balloon ports for ease of identifica Luer-lock balloon port to ensure secure syringe connection

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-14 to Fr-24	Fa	1	10	7453



ΔVΔNOS

MIC-KEY* Low-Profile Gastrostomy Feeding Tube

Avanos MIC-KEY* Low-Profile Gastrostomy Feeding Tube is designed to provide convenience to patients requiring gastrostomy feeding as it is easy to conceal, it sits on skin level, allowing limited interference to clothing.

- > Made from medical grade silicone
- > Low-profile design makes it unobtrusive and easy to conceal, sits at skin level
- Inflatable silicone, internal retention balloon reduces patient trauma
- > With included accessories for feeding and medication
- > With range of stoma length available

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	DVR
1-KC-LPGF12	Avanos MIC-KEY* Low-Profile Gastrostomy Feeding Tube, Fr-12	Ea	1	_	-
1-KC-LPGF14	Avanos MIC-KEY* Low-Profile Gastrostomy Feeding Tube, Fr-14	Ea	1	-	-
1-KC-LPGF16	Avanos MIC-KEY* Low-Profile Gastrostomy Feeding Tube, Fr-16	Ea	1	-	-
1-KC-LPGF18	Avanos MIC-KEY* Low-Profile Gastrostomy Feeding Tube, Fr-18	Ea	1	-	-
1-KC-LPGF20	Avanos MIC-KEY* Low-Profile Gastrostomy Feeding Tube, Fr-20	Ea	1	-	-
1-KC-LPGF24	Avanos MIC-KEY* Low-Profile Gastrostomy Feeding Tube, Fr-24	Ea	1	-	_





OrthoSenz[®] Orthopedic Cast Padding

Orthosenz Orthopedic Cast Padding is made from 100% cotton, designed to provide superior cushioning and efficient protection.

FEATURES AND BENEFITS:

- ► Made of non-woven, stretchable, 100% non-allergic cotton.
- ► Highly porous fabric assists in allowing circulation of air, reducing skin irritation.
- ► Easy to tear for conformability and easy application
- ► Durable to sustain wear and tear during use.

DESC	RIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available in: 3'	x 5 yards to 6"x 5 yards	Roll	10	400	00402




OrthoSenz[®] Plaster of Paris

Orthosenz Plaster of Paris has rich and creamy texture that provides excellent molding, smoother, stronger and appealing cast finish. Reliable cast that provides high comprehensive strength in short immersion time.

- Fast setting: 2-5minutes
- Lightweight, durable and adapts easily to the contours of the body.
- ► Non-toxic, non-allergic
- ► Individually packed in polybag

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: 4"x 5 yards and 6"x 5 yards	Roll	10	100	00401





Renasenz®

Blood Tubing Set

FEATURES AND BENEFITS:

- ► Medical grade materials displaying good biocompatibility
- Fluid pathway is smooth inside, preventing clotting or air bubbles
- ► Tubes are flexible, kink-resistant
- ► Color-coded tubings and ports
- ► Pump segment has excellent anti-fatigue performance
- ► Luer-lock connectors according to international standards, CE certified



REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-SM-BTS000	Renasenz Blood Tubing Set	Ea	1	25	09071

Renasenz®

Hollow Fiber Dialyzer RH18

High Flux

FEATURES AND BENEFITS:

- ► Micro-wave design ensures extended effective surface area
- ► Integrated Nano-spinning technology

- ► Unassigned vertical position: flexible to use dialyzer, place either which vertical direction.
- ► Exceptional clearance of Beta-microglobulin
- ► Housing Material: Polycarbonate, Membrane Material: Polysulfone
- ► Surface Area: 1.8m²

REF CODE	DESCRIPTION / SIZE	UOM	QTY/CTN	MDR
0-PM-DZHF18	Renasenz Hollow Fiber Dialyzer RH18	Ea	12	08712B



Renasenz®

Hollow Fiber Dialyzer R15, R18

Low Flux High Performance

FEATURES AND BENEFITS:

- ► Micro-wave design ensures extended effective surface area
- ► Integrated Nano-spinning technology
- ► Smooth cutting surface, reducing incidence of blood cell damage
- ► Smaller priming volume for cost-effective dialysis treatment
- ► Unassigned vertical position: flexible to use dialyzer, place either which vertical direction.
- ► Housing Material: Polycarbonate, Membrane Material: Polysulfone
- ► Surface Area: R18 (1.8m²), R15 (1.5m²)

REF CODE	DESCRIPTION / SIZE	UOM	QTY/CTN	MDR
0-PM-DZLF18	Renasenz Hollow Fiber Dialyzer R18	Ea	12	08712A
0-PM-DZLF15	Renasenz Hollow Fiber Dialyzer R15	Ea	12	08712A



Unimex®

AV Fistula with Back-Eye Needle

Unimex AV Fistula with Back-Eye Needle is made from high-quality material, designed to provide easy manipulation and effective dialysis treatment.

- ► Sharp-beveled needle minimizing patient trauma
- ► Flexible fixed wings, designed for perfect gripping during insertion
- ► Stainless steel needle with ultra thin wall
- ► Non-toxic, non-pyrogenic
- Available in G16 and G17

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available in G16 and G17	Ea	1	50	500	10111



Renasenz®

Hemodialysis Catheter Kit

Renasenz Hemodialysis Catheter kit contains a hemodialysis catheter that is made from superior quality Thermoplastic Polyurethane material with innovative catheter tip orientation, designed for optimal flow and easy maintenance during hemodialysis therapy.

FEATURES AND BENEFITS:

- ► Made from Medical Grade PU (Polyurethane), provides better flexibility and anti-kinking performance.
- Thermosensitive, softens according to the patient's body temperature.
 Smooth inside and outside surface to avoid clot formation.
- ► Radiopaque, facilitates quick visualization under X-Ray, for placement confirmation.
- Rotatable suture wings, facilitates skin inspection.
 Tapered tip facilitates ease of insertion and reduces vessel trauma during insertion.
- ► Silicone extension tube to increase patient comfort.
- Side holes are 2.5cm apart from the end hole to prevent wall-suction, minimizing the risk of clot formation.
 Inline clamps with printed priming volumes.
- ► Introducer needle with Pressure Activated Safety Valve
- Available in Double Lumen and Tri-Lumen: Standard and Compound Kit

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available in 2 and 3 Lumen, Curved and Straight	Set	1			10719



Unimex[®]

Blood Bag

Unimex Blood Bag is manufactured from high-quality medical grade PVC sheet with tamper evident label, available in single, double and triple blood bag system; designed for aseptic and convenient collection and separation of whole blood and blood components.

FEATURES AND BENEFITS:

- ► Siliconized ultra-thin wall G-16 sharp beveled needle assures smooth an atraumatic venipuncture
- ► Clear, soft, frosted and flexible kink-resistant tubing with unique identification number at intervals of
- ► Tamper evident label allows easy marking/writing that cannot be removed or tampered with.
- ► Double and Triple Blood Bag System are designed with snap tip and connector tube to allow smooth transfer of blood components and additive solutions.
- ► Convenient hanger slits allows easy suspension of the bags during collection and storage.
- ► Anticoagulant CPDA solution preserves whole blood up to 35 days

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in: Single, Double and Triple Bag	Ea	10	100	08471A



Unimex®

Blood Transfusion Set

Unimex Blood Transfusion Set is designed for convenient administration of blood and blood components.

- ► Sharp spike for better penetration of plastic bottle
- ► Dripping port provides constant precision drip over long duration
- ► 20 drops/ml for accurate flow rate
- ► Blood filter as per ISO standards
- ► Mesh filter is made from monofilament fiber for effective filtration
- ► With Y-injection port for easy administration of medication
- ► 200cm kink-resistant transparent tubing which provides clear view of potential blockage or backflow

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-UN-BT002C	Unimex Blood Transfusion Set, 200cm	Ea	20	400	07418





Respisenz®

Oxygen Connecting Tube

Respisenz Oxygen Connecting Tube is designed to provide easy connection and pathway of uninterrupted oxygen flow from the tank to the patient receiving oxygen therapy.

FEATURES AND BENEFITS:

- ► Made from soft lightweight medical grade PVC
- ► With distinctively designed 5 channel inner lumen for continuous airflow
- ► Crush-resistant
- ► 2.1 meters in length
- ► Latex-free





Respisenz[®]

Nasal Oxygen Cannula

Respisenz Nasal Oxygen Cannula is designed for long term oxygen treatment. It carries 1-5 LPM of oxygen and delivers 28-44% of oxygen concentration, depending on the oxygen flow and patient's breathing pattern.

FEATURES AND BENEFITS:

- Made from soft, lightweight medical grade PVC, safe to use and hypoallergenic
- ► Smooth rounded surface and edges to prevent local irritation of the nasal mucosa
- ► With 2 meter tubing
- ► Available in Neonate, Pedia and Adult

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in: Neonate, Pedia and Adult	Ea	1	100	05607



Respisenz®

High Flow Nasal Cannula

A non-invasive respiratory support, connected to the oxygen delivery system, designed to deliver high flow of blended air and oxygen through the nasal interface. This device is capable of delivering up to 100% humidified and heated oxygen at a flow rate of up to 60 liters per minute.

- ► Prongs are made from silastic material, designed with larger cavity to reduce nasal resistance, obstruction and to ensure the delivery of dynamic positive airway pressure.
- ► Nasal prongs are gentle on the nose: improves patient comfort, compliance and tolerance to therapy
- ► Patented adapter that swivels 360° to ensure kink resistance of the tubing.
- Flexible tubing with silicone material, especially designed to address water condensation.
- ► Connector with multi-diameter interface for connecting different sizes of tubing.
- ► With Fisher & Paykel third generation pipe joint for compatibility connection.

	REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
i	0-PM-HFNCS0	Available in: Small. Medium and Large	Ea	10	100	_



Respisenz[®]

Medium Concentration Oxygen Mask with Tubing

Respisenz Medium Concentration Oxygen Mask with tubing is used to deliver oxygen to a patient requiring moderate concentration of oxygen for a short period of time.

FEATURES AND BENEFITS:

- ► Uses 6-8 LPM and gives 60-80% oxygen concentration
- ► With adjustable metal nose clip
- Has smooth rounded surface, preventing pressure points to the patient's face
 Packed with 7 feet crush-resistant tubing

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in: Pedia and Adult	Ea	1	100	05370



Respisenz[®]

High-Concentration Oxygen Mask with Tubing

Respisenz High-Concentration Oxygen Mask with Tubing is an oxygen mask with tubing and reservoir bag that provides increased oxygen supply needed by the patient. The reservoir bag provides adequate oxygen supply for successful oxygen treatment.

FEATURES AND BENEFITS:

- ► Made from high-quality medical grade PVC
- ► Designed with turn-up rim to ensure comfortable fit and good seal
- ► With 7 feet kink-resistant, star lumen tubing that ensures patient mobility and continuous oxygen flow
- ► Tubing adaptor swivels to accompdate patient's position
- ► Mask has smooth molded surface to prevent pressure points and trauma to the patient's face

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in: Pedia and Adult	Ea	1	100	05370



Respisenz®

Multi-Vent Mask

Respisenz Multi-Vent Mask is an oxygen mask with an adjustable venturi device (diluter) that entrains and mixes specific amount of oxygen with room air to attain precise oxygen concentration.

FEATURES AND BENEFITS:

- ► Made from high-quality medical grade PVC, safe-to-use and hypoallergenic
- ► Air-entrainment / jet nozzle with adjustable settings
- ► 7 feet tubing with star lumen to allow continuous air-flow even when kinked
- ► Clear and soft mask with pliable nose clip to ensure comfortable fit while allowing observation of the patient's face
- ► Tubing adaptor swivels to accomodate patient's position

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in: Pedia and Adult	Ea	1	100	08320



Respisenz[®] **Aerosol Mask**

Respisenz Aerosol Mask is an elongated mask that fits under the chin, designed to provide patient comfort during aerosol therapy.

- ► Hypoallergenic elastic band with adjustable metal nose clip for better fitting
- ► With dual exhalation port
- ► Smoothly molded edges, prevents pressure and trauma to the patient's face
- ► Clear and transparent, allows observation of the patient's face

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in: Pedia and Adult	Ea	1	200	05370



Respisenz®

Nebulizer Kit with Mouthpiece

Respisenz Nebulizer Kit with Mouthpiece is designed for hand-held aerosol therapy. It delivers complete medication producing an exceptionally fine mist, consistent in particle size.

FEATURES AND BENEFITS:

- ► Made from high-quality polyethylene material, durable and safe to use
- With anti-spill chimney and packed with 7 feet crush-resistant tubing
 Operates at angles to accommodate common patient position while maintaining treatment
- ► Delivers high respirable fraction of greater lung deposition

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Respisenz Nebulizer Kit with Mouthpiece	Ea	1	50	05474



Respisenz®

Nebulizer Kit with Mask

Respisenz Nebulizer Kit with Mask is designed for in-line aerosol therapy. It delivers complete medication, producing an exceptionally fine mist, consistent in particle size.

FEATURES AND BENEFITS:

- ► Made from high-quality polyethylene material, durable and safe to use
- Delivers high respirable fraction of greater lung deposition
 Operates at angles to accomodate common patient position while maintaining treatment
- ► Anti-spill design to prevent medication waste even at high pressure
- ► With 7 feet crush-resistant tubing

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in: Pedia and Adult	Ea	1	100	05474



Respisenz®

Pre-Filled Humidifier with Adaptor

Respisenz Pre-filled Humidifier with Adaptor is used during oxygen therapy to provide humidification, removing dryness in oxygen, reducing mucosal damage, offering comfort and safety to the patient while preventing cross-contamination.

FEATURES AND BENEFITS:

- Filled with sterile water for inhalation, prepared by distillation
- Combines an easy to use adapter which pierces into the inlet sterile bottle to secure attachment
- ► Has multi-angle outlet port that prevents spitting even in high flow rates
- With small bottle opening that prevents refilling
 4 psi pressure-relief valve emits audible alarm when line is occluded or in high-pressure
- ► Bottle is made from high-density polyethylene

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-USH350	Respisenz Pre-Filled Humidifier with Adaptor, 350ml	Bottle	1	25	03999



Respisenz®

Pre-filled Humidifier 1000ml

Respisenz Pre-filled Humidifier 1000ml is designed for humidification during inhalation therapy, removing dryness in oxygen to reduce mucosal damage, offering comfort and safety to the patient.

- ► Made from sterile purified water
- ► Free from antimicrobial substances
- Allows simplified replenishing empty refilling chambers and nebulizers when attached to appropriate feed sets.
- ► One-step set-up helps reduce patient interruptions
- ► Sterile water in hanging bag with 1000ml volume
- ► Hanging bottles with sterile water for inhalation features a closed system design to help minimize contamination

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-HB1000	Respisenz Pre-Filled Humidifier, 1000ml	Bag	_	20	03999



Respisenz®

Manual PVC Resuscitator

Respisenz Manual PVC Resuscitator is a handheld device, used to manually ventilate apneic patients, also used for spontaneously breathing patients to augment ventilation and/or oxygen delivery.

FEATURES AND BENEFITS:

- ► Single-use to prevent cross-contamination
- ► Comes with clear mask for immediate use
- ► ISO standard connectors
- ► Fast recoil time for rapid ventilation
- ► Equipped with pressure relief valve for easy adjustment of pressure requirement
- ► Made from medical grade PVC
- ► Available in pre-assembled type to provide dependability during emergency situations

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available in: Infant, Pedia and Adult	Set	_	1	24	05294



Unimex®

Suction Catheter

Unimex Suction Catheter is used for oral and tracheo-bronchial suctioning of secretions to maintain airway patency.

FEATURES AND BENEFITS:

- Rounded distal eye with 2 lateral relieving eyes that reduces risk of blockage and trauma to the patient's tissue
- $\,\blacktriangleright\,$ Frosted medical grade PVC material allows atraumatic introduction
- ► Also available in suction catheter with one piece PE Glove

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-05 to Fr-18	Ea	1	50	500	05746




simplex

Suction Catheter

Simplex Suction Catheter is a thin sterile tube designed to draw out secretions such as mucus or saliva in the upper respiratory tract. With frosted surface that provides low friction, making it atraumatic to the patient's delicate mucosa.

- ► Smooth rounded tip and lateral eyes prevents trauma to the patient's tissue
- ► Color-coded for size identification
- ► Frosted surface allows smooth introduction of the catheter
- ► Non-toxic and kink-resistant

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-05 to Fr-18 (53cm)	Set	1	50	500	05320





Suction Catheter with Finger Tip Control

Unomedical Suction Catheter with Finger Tip Control is used to remove secretions from the mouth and airways of critically ill patients to maintain airway patency.

FEATURES AND BENEFITS:

- ► Color-coded for easy catheter identification
- ► Frosted surface facilitates easy and smooth introduction
- ► Transparent tubing for observation of secretions

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-05 to Fr-18	Ea	1	75	450	_



Respisenz[®]

Closed Suction Catheter System

Respisenz Closed Suction Catheter System allows safe suctioning of secretions and continuous ventilation during suctioning, designed to prevent VAP, PEEP loss, nosocomial infection and cross-contamination. With built-in MDI port that provides accessible aerosol medication administration.

FEATURES AND BENEFITS:

- ► 100% sealed connection to prevent cross-contamination and the spread of infection
- ► Double-swivel elbow connector to provide patient comfort
- ► Frosted catheter with atraumatic tip and lateral eyes, reducing mucosal damage
- ► MDI port for drug delivery and flushing port for cleaning of the catheter
- ► 72 hours maximum usage with day sticker for replacement guide

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available Size: Fr-10 to Fr-16	Set	1	20	00086



ΔVΔΝΟ5

Turbo Cleaning Closed Suction System (Endotracheal)

Avanos Turbo Cleaning Closed Suction System is intended to safely remove secretions from the trachea and bronchi while maintaining oxygen supply to a patient with endotracheal tube.

- ► Double swivel elbow for suctioning and ventilation
- ► Cleaning chamber with twin peep seal technology
- One-way lavage port prevents spray-back of saline during catheter cleaning
- ► 72 hours maximum use
- $\,\blacktriangleright\,$ With turbulent cleaning action, decreases bacterial colonization on the catheter tip by 89%

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available Size: Fr-12 to Fr-14	Ea	1	20	7454



ΔVΔNOS*

Trachcare Closed Suction Catheter System (Tracheostomy)

Avanos Trachcare Closed Suction Catheter System is intended to safely remove secretions from the trachea and bronchi while maintaining oxygen supply and ventilation to a patient with tracheostomy tube.

FEATURES AND BENEFITS:

- ► Double swivel elbow for suctioning and ventilation
- ► Cleaning chamber with twin peep seal technology
- $\buildrel {}^{\buildrel {}^{$
- ► Thumb valve with 360° locking system prevents accidental suctioning
- ► Turbulent cleaning action decreases bacterial colonization on the catheter tip by 89%

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available Size: Fr-12 to Fr-14	Ea	1	20	7454



©Unomedical

Cathy Closed Suction System

Unomedical Cathy Closed Suction System is a closed system catheter, designed to prevent cross-contamination and Ventilator Associated Pneumonia (VAP).

FEATURES AND BENEFITS:

- ► Unique disconnecting ring for easy disconnection, preventing tube dislodge
- Swivel standard connector that fits all standard ventilator equipment, allowing suctioning and ventilation at the same time
- ► Convenient knob that clearly indicates open or close position
- ► Clear safety lid that prevents accidental press of the suction control button.
- ► Color-coded flanges for immediate visual identification of sizes
- ► Universal connectors for attachment to all suction tubes and sputum sample collection devices

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available Size: Fr-12 to Fr-14	Ea	1	18	8027





Mucus Extractor

Simplex Mucus Extractor is a mucus specimen container with universal tubings and an extra lid, designed for obtaining mucus specimen for examination. Available in two designs: (1) Direct to Patient: with integrated suction tubing that goes directly to the patient, minimizing the use of additional suction catheter. (2) With Catheter Connector: designed with funnel connector that provides universal connection to any size of open/closed suction catheter.

- $\buildrel {}^{\buildrel {}^{$
- Extra lid ensures safe sealing of specimen container
- ► Tubings are made from soft and kink-resistant material
- ► Universal connectors for easy connection
- ► 25ml container with graduation markings for quick volume detection

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
2-PP-IME001	Simplex Mucus Extractor DP w/ 50cm Tubing, Fr-10	Ea	1	25	200	_
2-PP-IME002	Simplex Mucus Extractor w/ Catheter Connector	Ea	1	25	200	_



Pulmodyne[®]

Soft Touch Tracheostomy Tube Holder

Pulmodyne Soft Touch Tracheostomy Tube Holder is especially designed to provide comfortable support, while securing tracheostomy tube positioning, thus preventing accidental dislodgement and movement of the tube.

FEATURES AND BENEFITS:

- $\begin{tabular}{ll} \begin{tabular}{ll} \beg$
- Latex-free, foam padded with soft cotton lining
- Provides snug fit to tracheostomy tube while providing patient comfort

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-BL-313720	Pulmodyne Soft Touch Tracheostomy Tube Holder	Ea	20	200	MDRE



Respisenz®

Tracheostomy Mask

Respisenz Tracheostomy Mask is designed for Laryngectomy or Tracheostomy aerosol therapy.

FEATURES AND BENEFITS:

- ► Tubing connector swivels 360° for ease of positioning
- ► 22mm OD connector for corrugated tubing and nebulizer tees
- Exceptional strap clip allows easy mask detachment to lessen patient discomfort

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-SM-TM001A	Respisenz Tracheostomy Mask, Adult	Ea	1	50	05370



Panamed®

Maintenance Care Pack

Panamed Maintenance Care Pack offers the finest in nursing convenience and infection control, this pack provides all the necessary components, needed for efficient and quick tracheostomy/endotracheal care procedures.

FEATURES AND BENEFITS:

- ► Assembled in a compact kit
- ► Peel-down lid protects contents during storage and handling
- ► Efficient, convenient to use and cost economical

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-PP-TRCKIT	PANAMED Maintenance Care Pack	Kit	1	50	_



Respisenz®

Spirometer Filter with Nose Clip

Respisenz Spirometer Filter with Nose Clip is used during pulmonary function test, designed to provide barrier filtration for protection of patient, staff and equipment from bacterial and viral

- ► Filtration: Electrostatic / Mechanical
- ► Flow Resistance: 0.93 cmH20
- ► Bacterial Filtration Efficiency: >99.99%
- ► Viral Filtration Efficiency: >99.99%
- ► Filtration Area: 11.46cm
- ► Exchange Filtration Area: >34.00cm
- ► Oval Internal Dimension: 18.00 x 34.00mm
- ► Weight: 18g

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-SM-SPFTNC	Respisenz Spirometer Filter with Nose Clip	Kit	1	90	_



Respisenz®

Heat and Moisture Exchanger (HME) Filter

Respisenz Heat and Moisture Exchanger (HME) Filter is a type of filter with humidification properties, aside from it's filtering capacity, used with patients under mechanical ventilation, designed to prevent drying of the respiratory mucosa.

FEATURES AND BENEFITS:

- ► Filtration Efficiency (BVF & VFE): 99.99%
- ► Humidity output: 31mg/H20 @VT500ml
- ► Resistance: 180pa @30LPM
- ► Dead Space: 45ml ► Tidal Volume: 150-1500ml
- ► Connections: 22M / 15F-22F / 15M
- ► Weight: 30g
- ► Gas monitoring port: (luer port) with retainer strap

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-SM-BVF001	Respisenz Heat and Moisture Exchanger (HME) Filter	Ea	1	100	_



Respisenz[®]

Bacterial Viral Filter (BVF)

Respisenz Bacterial Viral Filter (BVF) is generally consisting of a flat wad of electrostatically charged fibers, designed to trap bacteria and viruses in a closed breathing system.

FEATURES AND BENEFITS:

- ► Bacterial Filtration Efficiency: 99.99%
- ► Viral Filtration Efficiency: 99.99%
- ► Resistance: 60pa @30 LMP
- ► Dead Space: 32ml
- ► Tidal Volume Range: 250-1500ml
- ► Gas monitoring luer port with retainer strap
- ► Connections: 22M/ 15F-22F/ 15M

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-SM-BVF002	Respisenz Bacterial Viral Filter (BVF)	Ea	1	100	_



Respisenz®

MDI Adapter

Respisenz MDI Adapter is a connector that allows delivery of metered dose inhaler to a patient hooked in a mechanical ventilator.

FEATURES AND BENEFITS:

- ► Made from Polyethylene material
- Accepts plastic or metal tipped standard MDI canister
- Arrow marker for the direction of flow
- Available in straight and T-connector
- ► OD: 22mm. ID: 22mm

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-MDI001	Respisenz Metered Dose Inhaler (MDI) Connector, Straight	Ea	1	500	889438
0-PM-MDI002	Respisenz Metered Dose Inhaler (MDI) T-connector	Fa	1	500	889438



Respisenz®

Corrugated Aerosol Tubing

Respisenz Corrugated Aerosol Tubing is a flexible tubing made from polyethylene and ethyl vinyl acetate plastic, designed for efficient aerosol therapy.

- ► Continuous 100ft roll of flexible, kink-resistant tubing
- ► Segmented with cuttable sections every 6 inches
- ► Collared tubing end slips, allowing secure connection
- ► Clear tubing with cuff that allows easy universal 22mm connection

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-SM-CT0001	Respisenz Corrugated Aerosol Tubing, 100 Feet/Roll	Roll	1	1	_



Respisenz®

Ventilator Circuit

Respisenz Ventilator Circuit is an assembly of components which connects the patient's airway to the mechanical ventilator, providing oxygen to a patient who is physically unable to breath or is breathing insufficiently.

FEATURES AND BENEFITS:

- ► Bicolored corrugated tubings, allowing identification of either the inspiratory or the expiratory limb
- ► Standard connectors to fit with most ventilators and humidifiers
- ► Flexible and lightweight
- ► With available options to accomodate different patient requirement
- ► 1.8 meter breathing circuit adult manifold

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in Nebulizer, Water Trap and Universal (Adult)	Set	1	20	



Respisenz®

Peak Flow Meter

Respisenz Peak Flow Meter is a portable handheld respiratory device that measures a person's ability to push air out of the patient's lungs. Airflow is measured by the amount of air that can blow out in one "fast blast". It helps in managing and monitoring the patient's asthma by detecting airway narrowing.

FEATURES AND BENEFITS:

- ► Mouthpiece size is convenient to use for pedia and adult, with easy-to-grip handle
- ► Clearly labeled numerical scale with indicator arrows to easily monitor asthma control
- ► With monitoring table for tracking

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-UP1038	Respisenz Peak Flow Meter	Ea	1	50	841068



Respisenz[®]

Incentive Respiratory Exerciser

Respisenz Incentive Respiratory Exerciser helps patient to recover normal respiration after a chest or abdominal surgery by encouraging deep and prolonged inspiration.

- Designed with three chambers to provide patient with smooth breathing exercise stepwise from 600ml, 900ml and 1200ml
- ► Compact design and cost economical
- ► Single-patient use

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
3-GC-MIS000	Respisenz Incentive Respiratory Exerciser	Ea	1	28	773866



Lifotronic

High-Flow Humidification System

A humidifier with an integrated flow generator that delivers HIGH-FLOW warmed and humidified respiratory gas to spontaneously breathing patients through a variety of interface. Designed to deliver high FiO2, to accurately meet the patient's ventilatory demands, while providing patient comfort.

FEATURES AND BENEFITS:

- ► Delivers a modest amount of positive airway pressure.
- ► 5-inch colored touchscreen with simple user interface
- Multiple temperature settings (31-37°C, increment 1°C) to improve comfort and compliance of treatment.
- ► Built-in oxygen blender. Oxygen concentration can be adjusted from 21 to 100%.
- ► One-way valve design to avoid contamination of the unit.
- ► With twelve alarm functions
- ► With available replacement filter and disinfection tube to carry out disinfection cycle



PARTS OF THE SYSTEM:



Heated Humidifier



Disposable Patient Circuit Tube and Chamber Kit






High Flow Nasal Cannula Flow range: Small 10-50L/min Medium/Large 10-60L/min



Air Filter



REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-LT-HFH660	Lifotronic High-Flow Humidification System	Unit	_	1	832319



WSTERIPAK®

Bowie and Dick Sterilizer Efficiency Test Pack

Steripak Bowie and Dick Test Pack is a chemical indicator which contains a series of air removal and steam penetration barriers. Successful test confirms that steam penetration into the test pack is rapid and even.

FEATURES AND BENEFITS:

- ► Prevents problems associated with imperfect heat sterilization
- ► Helps identify sterilization leakage and efficiency
- ► Chemical indicator inside the pack color changes from yellow to black
- Non-toxic indicator sheet, can be kept as a permanent record
- ▶ Lead-free

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
2-PP-BODKTR	Steripak Bowie and Dick Test Pack	Pack	25	4	100	_



WASTERIPAK®

Steam Autoclave Sterilization Indicator Tape

Steripak Steam Autoclave Tape is a specially formulated sterilization adhesive tape which serves as a reliable chemical indicator between processed and unprocessed item.

FEATURES AND BENEFITS:

- Adheres strongly to woven or non-woven paper and plastic wraps
- ► Securely adheres under high temperature and high pressure conditions
- ► Indicator lines change without ink rub off

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available in 3/4", 1/2"and 1"	Roll	1			_



Unimex®

Autoclave Tape

Unimex Autoclave Tape is a especially formulated sterilization indicator tape, which serves as a reliable chemical indicator whether the pack is processed or not.

FEATURES AND BENEFITS:

- ► Backing made from saturated crepe paper
- ► Natural rubber, non-staining, adhesive spread
- ► Indicating lines change from white to black without ink rub off
- ► Adheres to woven or non-woven, paper and plastic wraps

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available in 3/4", 1/2"and 1"	Roll	1			_



WSTERIPAK®

Self-Sealing Sterilization Pouch

Steripak Pre-Cut Self-Sealing Sterilization Pouch is available in various sizes that caters low profile to medium instruments devices. Constructed from high-quality medical-grade film and paper suited for sterilization. Provides quick and secure closing of packs without a heat sealer.

- ► Made from 60gsm medical grade paper for superior barrier protection against contamination
- ► Bioadhesive strip that ensures seal integrity
- $\begin{tabular}{ll} \begin{tabular}{ll} \beg$

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
2-PP-SPSS57	57x130mm / 2.24x5.11"	Box	200	30	_
2-PP-SPSS70	70x260mm / 2.75x10.23"	Вох	200	15	_
2-PP-SPSS90	90x260mm / 3.54x10.23"	Box	200	10	_
2-PP-SPS135	135x280mm / 5.31x11.02"	Вох	200	10	_
2-PP-SPS190	190x360mm / 7.48x14.17"	Box	200	5	_



WSTERIPAK®

Sterilization Gusseted Reel


Steripak Sterilization Gusseted Reels are made from medical grade paper and multi-layer barrier film, designed for easy and reliable sterilization. Ideal for all packaging demands of medical instruments and kits

FEATURES AND BENEFITS:

- ► Superior barrier with 60gsm medical grade paper
- ► Gusset folded on each side of the pouch which enables a bigger inner space than the flat pouches
- ► Minimized printed area in and between the sealed layer of paper and film to prevent contamination

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in 3" to 16" (100m)	Roll	1	8	_





WSTERIPAK®

Sterilization Flat Reel

Steripak Sterilization Flat Reels are made from medical grade paper and multi-layer barrier film, designed for easy and reliable sterilization. Ideal for small and lightweight items/devices.

FEATURES AND BENEFITS:

- ► Made from 60gsm medical grade paper
- ► Flat type, perfectly fits to all lightweight packaging
- ► Reinforced film provides tear resistant quality
- ► Printed with double process indicators for steam and EO sterilization

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
0-PM-LPEG1L	Steripak Fat Reel, 50mm x 200m (2") x 656ft		Roll	1	12	_



WSTERIPAK®

Tyvek® Sterilization Flat Reel

Steripak Tyvek $^{\circledR}$ Sterilization Flat Reels are made from medical-grade paper and multi-layer barrier film, designed for easy and reliable plasma sterilization.

- Made from 60gsm medical grade paper, sealed on 2 sides with non-woven structure of flash spun
 polyethylene fiber
- ► High-density polyethylene, providing excellent strength and peelability
- ► Tear-resistant for superior microbial barrier
- ► Suitable for plasma sterilization

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-PP-ST1520	Steripak Tyvek $^{ ext{(B)}}$ Flat Reel, 150mm x 200m (6") x 656ft	Roll	1	4	_



WASTERIPAK®

Steam Indicator Strip

Steripak Steam Indicator Strips are chemical indicator strips that provides distinct color change when exposed to steam sterilization process. Placed in every pack, pouch, wrap and tray to distinguish sterilization failures through an immediate result.

FEATURES AND BENEFITS:

- ► Lead-free ink
- $\,\blacktriangleright\,$ Strips can be written on and kept as a permanent record
- Provides immediate identification (completion is indicated when the entire length of the strip turns from yellow to black when exposed to steam sterilization)
- ► Complies with EN-ISO 11140-1
- ► Class 4 indicator strips

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-ST-TSAI25	Steripak Steam Indicator Strip	Box	250	80	_



STERIPAK®

Ethylene Oxide (EO) Indicator Strip

Steripak Ethylene Oxide (EO) Indicator Strips provides distinct color change when exposed to EO sterilization process. Placed in pack, pouch, wrap and tray to distinguish sterilization failures through an immediate result.

FEATURES AND BENEFITS:

- ► Precise response to critical parameters of EO sterilization
- ► Indicator turns from red to blue when exposed to EO sterilization
- ► Complies with EN-ISO 11140-1
- ► Strips can be written on and kept as a permanent record
- ► Provides immediate identification
- > Class 4 indicator strips

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-PP-TSE025	Steripak Ethylene Oxide (EO) Indicator Strip	Box	250	80	_



WSTERIPAK®

Plasma Indicator Strip

Steripak Plasma Indicator Strips provides distinct color change when exposed to Plasma sterilization procedure. Placed in pack, pouch, wrap and tray to distinguish sterilization failures through an immediate result.

- ${\color{red} \blacktriangleright} \ \, \text{Internal indicator for all packages sterilized in Gas Plasma Hydrogen Peroxide sterilization system}$
- ► Color changes from green to yellow after successful sterilization
- ► Complies with EN-ISO 11140-1
- ► Strips can be written on and kept as a permanent record
- ► Provides immediate identification
- ► Class 4 indicator strips

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-PP-TSPL25	Steripak Plasma Indicator Strip	Box	250	80	_



W.STERIPAK®

Crepe Paper Wrapping Sheet

Steripak Crepe Paper Wrapping Sheets are effective bacterial barrier, used for inner or outer sterilization wrapping for various sets of medical devices, suitable for all types of sterilization methods to prevent infection and contamination.

FEATURES AND BENEFITS:

- ► Soft special paper, made from bleached, virgin wood pulp
- ► No odor, fiber-free and non-toxic
- ► Suitable for high-temperature, high-pressure plasma and ethylene oxide sterilization methods
- ${\color{red} \blacktriangleright}$ Made of 60gsm medical grade crepe paper material

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-PP-CP606	Steripak Crepe Paper, 60cm x 60cm, Blue	Pack	1	500	_



STERIPAK®

SMS Wrapping Sheet

Steripak SMS Wrapping Sheets are effective bacterial barrier, used as an inner or outer sterilization wrapping for various sets of medical devices, suitable for all types of sterilization methods to prevent infection and contamination.

- The SMS material is combined with multi-layer spunbond, meltblown and spunbond layer of polypropylene
- ► Extremely durable material with high tear resistance and burst strength
- ► Use as such or in combination with crepe papers
- ► Made of 60gsm medical grade paper
- ► Offers good drapeability, eliminates springback of sheets

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-PP-SWS120	Steripak SMS Wrapping Sheet, 120cm x 120cm	Pack	1	200	_





Panamed®

Disposable Surgical Gown

Panamed Disposable Surgical Gown are high-quality disposable fabric, designed to provide protection against infection during medical and surgical procedures. Offers superior protection from fluid exposure during invasive and intensive-manipulative procedures.

FEATURES AND BENEFITS:

Made from fluid-repellent, non-woven fabric with micro-fiber that provides high-resistance from fluid and bacteria

Abrasion-resistant and low-lint fabric to reduce risk of wound contamination post-operatively

€otton white knitted cuff sleeves for better fit and easy donning of gloves

∀elcro on the neck and ties at the waist for secure wearing

∀arious sizes available for superior fit and flexibility

Available in Basic, Standard and Premium

With 2 hand towels

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in Basic: Medium, Large and XL	Pack	1	30	752464
Available in Standard: Medium, Large and XL	Pack	1	30	752464
Available in Premium: Medium, Large and XL	Pack	1	30	752464





Surgical Gown

Vi Pha Co Surgical Gowns are specifically designed to provide efficient protective barrier against fluid and microbial transmission during medical and surgical procedures. Available in two types: Standard and Standard Reinforced.

- Made from multiple layers, providing efficient resistance against fluid and microorganisms, while providing high level of comfort to the wearer
- ► Constructed from abrasion resistant, air permeable, anti-static and low-linting fabric
- The Standard Reinforced is designed with specialized fabric reinforcement on chest and sleeves for maximum fluid protection
- ► Packed with two hand towels

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in Medium, large nd XL	Ea	1	40	00240



Standard

Standard Reinforced

Panamed®

Bouffant Doctor's Cap

Panamed Bouffant Doctor's Cap is an air-permeable, non-woven surgical cap which offers scalp and hair protection, with ties for secure wearing minimizing risk of falling hair into the sterile area during surgery.

FEATURES AND BENEFITS:

- ► Constructed with soft, durable ties for secure fitting
- ► Has a reinforced stitching for extra strength
- ► One size fits all

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
2-PP-DC25B0	Panamed Bouffant Doctor's Cap, Blue	Box	100	10	MDRE



Panamed®

Bouffant Nurse Cap

Panamed Bouffant Nurse Cap are sanitary head covers, providing protection for cleaner work environment. The breathable barrier design makes it suitable for hot and humid environment.

FEATURES AND BENEFITS:

- ► Elastic headband offers a comfortable, non-slipping fit
- Lightweight breathable fabric
- ► One size fits all

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available in Size 21" and 24"	Box	100	20	MDRE



Panamed®

Shoe Cover (Non-Skid)

Panamed Shoe Cover, Non-Skid is designed for comfortable wear in different sterile areas of the hospital like the operating room and delivery room.

FEATURES AND BENEFITS:

- ► Made from non-woven, lightweight fabric
- ► Waterproof and dust proof
- $\ \ \ \ \$ Anti-Skid treads on the bottom for safety

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-PP-OSCB00	Panamed Shoe Cover, Non-Skid, Blue	Pack	100	10	MDRE



Unimex®

Sterile Surgical Gloves

Unimex Sterile Surgical Gloves is designed to provide hand protection for the benefit of the patient and healthcare professionals alike during surgical operations.

- ► Powder-free and Hypoallergenic
- ► Tough glove yet gentle on hands
- ► Thin enough to provide sensitivity, thick enough to protect hands from contaminating materials

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available Size: 6.0 to 8.0	Pair	1	50	500	07870



Unimex[®]

Angiography Drape

Unimex Angiography Drape is a sterile drape with femoral and radial holes, designed for convenient use during Angiography procedures.

FEATURES AND BENEFITS:

Extended length and improved absorbent fabric reinforcement to minimize risk of fluid pooling


tow-lint fabric reinforcement to help reduce the risk of lint particles

One piece radial/femoral drape reduces waste and eliminates the use of additional drape

Radial entry advantage includes less bleeding to reduce risk for contamination

Available in two designs

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-600016	Unimex Angiography Drape (Lulu)	Ea	1	5	_
0-PM-600018	Unimex Angiography Drape (Orlando)	Ea	1	5	_



Unimex[®]

Laparotomy Pack

Unimex Laparotomy Pack is a complete procedural pack, designed to create a protective barrier between sterile and non-sterile field during laparotomy procedures.

FEATURES AND BENEFITS:

Made from Polypropylene SMS drape fabric, soft and drapeable while providing superior barrier strength

Fluid-resistant and impervious to minimize the risk of surgical site infection

Absorbent reinforcement fenestration surrounding the surgical site area to prevent fluid pooling during surgery

Full body coverage ensures aseptic technique

Saves time in pre-op collection of materials and unwrapping of single use products

Provides medical staff protection from bloodborne pathogens and other potential infectious agent

With five pieces breathable and fluid repellent surgical gowns

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-UVG00L	Unimex Laparotomy Pack	Pack	1	5	00232





Laparotomy Pack

VPC Laparotomy Pack is a complete sterile procedural pack, specifically designed with clinician's preference to obtain efficient barrier protection and infection control during laparotomy procedures.

FEATURES AND BENEFITS:

Made from impervious material, which prevents strikethrough of fluids, reducing risk of wound contamination and the associated risk of surgical site infection.

Aperture drape with optimal fluid control and reinforcement to ensure reliability and protection.

Additional utility drape for extra absorption.

Provides high barrier performance.

Puncture, tear and abrasion resistant.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-VP-PC0003	VPC Laparotomy Pack	Pack	1	3	00239







Obstetrics (OB) Pack

VPC Obstetrics (OB) Pack is a sterile pack that contains all the necessary components, essential during normal spontaneous delivery, designed with high level fluid collection for efficient infection control, for the well-being of both mother and child.

- Absorbent layer and impervious film for exceptionally high barrier strength and optimal fluid control.
- ► Drapes that are designed to sustain heavy fluid operation, with super absorbent reinforcement material in critical zone areas.
- ► Delivery pack with fluid collection
- Absorbent fabric that provides high strength and excellent drapability.
- Legging incorporated








Cesarean (CS) Pack

VPC Cesarean (CS) Pack is a sterile pack that contains all the necessary components essential during cesarean section delivery, designed with high level fluid collection for efficient infection control for the well-being of both mother and child.

FEATURES AND BENEFITS:

- Absorbent layer and impervious film for exceptionally high barrier strength and optimal fluid control.
- Drapes that are designed to sustain heavy fluid operation, with super absorbent reinforcement material in critical zone areas.
- Delivery pack with fluid collection.
- ► Absorbent fabric that provides high strength and excellent drapability.
- ► Highly efficient tape and incise film for longer operation. Small aperture for efficient infection control.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-PC0002	VPC Cesarean (CS) Pack	Pack	1	4	00164





Angiography Pack

VPC Angiography Pack is a sterile surgical pack, designed to provide efficient barrier protection and infection control during angiography procedures.

FEATURES AND BENEFITS:

- ► Drapes with optimal fluid control to ensure reliability and protection.
- ► Incise film provides effective barrier against microorganisms and bacteria, minimizing risk of infection.
- $\,\blacktriangleright\,$ Resistant to tearing and strikethrough.
- ► Reinforced absorbent material for optimal fluid protection.
- Four holes to access both radial and femoral arteries
- ► Transparent side panels to enable better visibility.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-VP-PC0006	VPC Angiography Pack	Pack	1	2	00062





Universal Pack

VPC Universal Pack is a sterile surgical pack, designed to provide simplified, reliable protection for many common surgical applications.

FEATURES AND BENEFITS:

- Composed of aseptic and high-quality components that provide safety, efficient fluid control, breathability and overall drape stretching resistance.
- Drapes with optimal fluid control to ensure reliability and protection.
- ► Super absorbent reinforcement material in critical zones for optimal fluid protection.
- Highly efficient tapes for long operations.
- ► Simple, flexible and adjustable drapes, adaptable to different common surgical procedures.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-VP-PC0004	VPC Universal Pack	Pack	1	7	00031





Ophthalmic Pack

VPC Ophthalmic Pack is a sterile surgical pack that offers efficient infection control, enhanced fluid protection and increased patient comfort in an ophthalmic or ocular surgery.

- Drapes with optimal fluid control, designed to effectively improve surgical efficiency and reduce the risk of cross-infection
- ► Incise film specifically designed for long operation
- Lightweight and breathable drape that provide efficient barrier protection, preventing fluid and microbial penetration.
- ► With fluid-collection pouch

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-VP-PC0005	VPC Ophthalmic Pack	Pack	1	40	00238



Unimex®

Gauze Ball with X-Ray Line (Peanut Sponge)

Unimex Gauze Ball with X-ray Line (Peanut Sponge) are gauze assembled in rounded shape, with a single X-ray element, designed to provide high absorbency and blunt dissection of delicate tissues.

FEATURES AND BENEFITS:

- ► Made from 100% cotton material
- ► With X-ray element
- ► 5cmx5cm, 28x24mesh, 4-ply

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
0-PM-BSX051	Unimex Gauze Ball w/ X-Ray Line (Peanut Sponge), 5cm x 5cm	Box	10	25	50	813693



Unimex[®]

Gauze Ball with X-Ray Line (Cherry Sponge)

Unimex Gauze Ball with X-ray Line (Cherry Sponge) are round tightly wrapped gauze with X-ray element, designed to provide high absorbency and blunt dissection of delicate tissues.

FEATURES AND BENEFITS:

- ► Made from 100% pure cotton material
- ► With X-ray line
- ► 10cmx10cm, 28x24mesh, 12-ply

REF CODI	E DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
0-PM-BSX101	Unimex Gauze Ball w/ X-Ray Line (Cherry Sponge), 10cm x 10cm	Box	5	20	20	813693



Unimex®

Gauze Swab/Sponge (Sterile & Non-Sterile)

Unimex Gauze Swab/Sponge products are made from 100% pure surgical grade cotton. Manufactured using a special process which produces softness and bulk for maximum absorption. Available in different sizes for wide range of dressing, cleaning, packing, prepping and debriding wounds.

FEATURES AND BENEFITS:

- ► Made from 100% pure surgical grade cotton gauze
- ► Soft and high in absorbency
- ► Folded edges and lint free
- ► Available with X-ray detectable line and no X-ray line
- > Pre-cut sizes ideal for user convenience

DESCRIPTION / SIZE Available Size in Sterile, Non-Sterile and with X-Ray Line: 4x4 and 4x8 (Mesh 28x24)



Unimex®

Absorbent Surgical Gauze Roll

Unimex Absorbent Surgical Gauze Roll are made from 100% cotton, ideal for making all types of sponges or dressing for different applications. Used for cleaning, packing, scrubbing, prepping and securing variety of wounds.

- ► Made from high-quality, highly absorbent natural pure cotton.
- Dyed in medical ways to ensure quality.
- ► Smooth, tucked in edges and no loose threads, excellent quality gauze.

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available Size: 20x12, 24x20 and 28x24	Roll	1	20	00022



Unimex[®]

Laparotomy Sponge

Unimex Laparotomy Sponge is an absorbent cotton fabric, used during surgical procedure. Provides excellent swab for bleeding control and tissue protection during surgical operations.

FEATURES AND BENEFITS:

- ► Made from 100% pure surgical grade cotton
- ► Manufactured using a special process which produces softness and bulk for maximum absorption
- High in absorbency and lint free
 All edges are folded in and stitched with cotton thread, preventing loose fibers

 The stiff of the st
- ► Radiopaque x-ray handle that remain visible even when saturated with blood

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in 18x18: 4 and 8ply	Pack	2	80	06710



Surgisenz®

Absorbable Gelatin Sponge

Surgisenz Absorbable Gelatin Sponge is a quick and effective hemostatic substance that stops bleeding fast, designed to reduce intra-operative bleeding and surgery time.

FEATURES AND BENEFITS:

- ► Allows uniform porosity which guarantees favorable hemostasis.
- Absorbs 40-50 times its own weight of water/oxalate blood.
 Completely absorbed within 3-4 weeks with no residue and no encapsulation, when implanted in vivo
- No cholesterol, pure gelatin, no preservatives, BSE free.
- ► Provides safety and excellent tolerability, no reactions or interactions with drugs or other food ingredients.
- ► pH neutral and Gamma-sterilized

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available in: 70x50 and 80x50	Ea	10	600	08257



Surgisenz®

Bone Wax

Surgisenz Bone Wax is a local hemostasis agent for the bone tissues; it is indicated for use in bleeding control by acting as a mechanical barrier during bone surgeries.

FEATURES AND BENEFITS:

- ► Sterile mixture of semisynthetic beeswax, paraffin and isopropyl palmitate
- ► Small amount is enough to achieve hemostasis
- ► Designed to minimize the delay of bone regeneration, possibility of foreign body formation and
- inflammatory reaction.

 Supplied sterile in 2.5 gram foil envelope with sealed overwraps.
- Sterilized by Gamma Radiation

REF CODE	DESCRIPTION / SIZE				MDR
0-PM-BW0025	Surgisenz Bone Wax, 2.5g	Box	12	100	10653



Panamed® Sterile Surgical Blade

Panamed Sterile Surgical Blade are used for procedures that requires puncture or cutting, for delicate dissection and critical tissue operation.

- ► Made from stainless steel material
- ► Known for precision and reliability
- ► Individually sealed in foil-peel sterile packaging for safety and ease of use
- ► Sterilized by Gamma Radiation (2.5m Rads)
- ► Available in multiple sizes
- ► 100 pieces per box

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: SS10 to SS24	Box	100	50	04837



Surgisenz®

Electrosurgical Cautery Pencil

Surgisenz Disposable Electrosurgical Cautery Pencil is used during surgical operations to make precise cuts and coagulate large bleeding areas.

FEATURES AND BENEFITS:

- Hexagon lock provides precision and reliability.
- ► Made from reliable stainless steel blade.
- ► Three pins banana plug, fits most electrosurgical generators.
- ► Length of ABS handle:160mm
- ► Length of Stainless Steel Blade: 48mm
- ► Length of Cable: 3.0m
- ► Diameter of Blade: 2.3mm

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-SM-ECP00	Surgisenz Electrosurgical Cautery Pencil	Ea	25	100	08455



Surgisenz[®] Electrosurgical Pad

Surgisenz Electrosurgical Pad is adhered to the patient's skin away from the surgical site, it is intended to provide a safe return path for electrosurgical current.

FEATURES AND BENEFITS:

- ► Single-use, avoiding cross-infection during surgery.
- ► Standard plug that works with most electrosurgical generators.
- ► Compatible with most skin types.
- ► Allows easy set-up
- Provides good drapability

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-SM-ESPB00	Surgisenz Electrosurgical Pad, Bipolar Adult	Ea	25	100	08403B



Surgisenz®

Cautery Tip Cleaner

Surgisenz Cautery Tip Cleaner is designed to effectively clean the tip of electrosurgical pencils and other surgical instruments.

FEATURES AND BENEFITS:

- Adhesive back for convenient placement during surgical procedures.
- ► Provides effective cleaning properties, does not cause damage to the tips of electrosurgical pencils.
- ► Non-abrasive material
- ► Provides high efficiency and durability
- ► Dimension: 5cm X 5cm

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-SM-CTC000	Surgisenz Cautery Tip Cleaner	Ea	100	400	_



Unimex®

Pre-Operational Disposable Razor

Unimex Pre-Operational Disposable Razor is used for shaving any length and density of hair before surgical operation.

- ► Blade made up of stainless steel material
- Well balanced plastic handle for optimal comfort
- ► Comes in a dispenser box

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-PM-UDRSS1	Unimex Single-Edge Disposable Razor	Box	100	10	MDRE



Unimex®

Skin Stapler - W35

Unimex Skin Stapler is used for skin closure as an alternative to sutures and needles. This device is easy-to-use and saves considerable time for the surgeon.

FEATURES AND BENEFITS:

- ► Angled head provides clear view for precise staple placement
- ► Staple count indicator clearly shows remaining staples
- Lightweight to reduce hand fatigue
- ► Excellent quality for reliable performance
- ► Shortens surgical time

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
1-UN-SS3500	Unimex Skin Stapler W35	Ea	1	10	80	04239



Unimex®

Disposable Skin Staple Remover

Unimex Disposable Skin Staple Remover is used for safe and quick removal of staples after wound healing.

FEATURES AND BENEFITS:

- ► Removes atraumatically all brands of surgical skin staples
- ► Provides fast and easy removal of staples
- ► Easy-to-handle making staple removal simpler

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
1-UN-SRR000	Unimex Disposable Skin Staple Remover	Ea	1	10	300	05727



Unimex®

Umbilical Cord Clamp

Unimex Umbilical Cord Clamp is used for clamping of umbilical cord of a newborn baby immediately after birth.

FEATURES AND BENEFITS:

- ► Provided with finger grip for safe and convenient handling
- ► Grooved clamping area to prevent slipping of the umbilical cord
- ► Secure lock to prevent accidental opening after clamping

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
0-PM-UUCC00	Unimex Umbilical Cord Clamp	Box	1	100	20	05726



☼ Unomedical

Umbilical Cord Clamp Clipper

Unomedical Umbilical Cord Clamp Clipper is a disposable device, used for the removal of an umbilical cord clamp.

- ► Disposable, non-sterile
- ► Compatible with most clamps
- ► Cuts hinge to remove umbilical cord clamp

REI	F CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
1-UM-0	06182 Unom	edical Umbilical Cord Clamp Clipper	Ea	1	25	200	_





Frazier Tip with Suction Connecting Tube

Simplex Frazier Tip with Suction Connecting Tube is designed for precise and extreme removal of fluids and debris from confined surgical sites.

FEATURES AND BENEFITS:

- ► Malleable shaft allows precise angle adjustment for accurate suctioning
- With stylet for easy declogging
- ► Finger control vent on the suction handle for better suction control
- ► With 210cm universal and kink-resistant tubing

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in Fr-8 to Fr-12	Ea	1	50	05111



simplex

Yankauer Suction Set

Simplex Yankauer Suction Set is used during surgical procedures to clear the operative site from fluid, blood and debris. Commonly used for suctioning small areas such as throat, arms or mouth.

FEATURES AND BENEFITS:

- ► Transparent tube allows clear visualization of fluids
- ► Sturdy, shatter-resistant construction allows effective suctioning
- ► Comes with 210cm kink-free transparent tubing
- ► With atraumatic relieving eyes, reducing the chances of tissue grab

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-SSYKRS	Simplex Yankauer Suction Set	Ea	1	50	01661





Yankauer Suction Set

Unomedical Yankauer Suction Set is intended for the removal of blood and fluids during surgery.

FEATURES AND BENEFITS:

- ► Made from high-quality, medical grade material
- ► Atraumatic tip minimizing tissue damage during suctioning
- ► 210cm non-collapsible tubing even at higher flow rate
- ► Funnel connector ensures tight grip

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-UMSY303	Unomedical Yankauer Suction Set 303	Ea	1	50	04258



simplex

Abdominal Poole Drain Suction Set

Simplex Abdominal Poole Drain Suction Set provides fast, efficient and safe suctioning of fluid/blood during surgical operations.

- ► Made from medical grade PVC
- ► Atraumatic distal tip minimizing tissue damage
- ► Tip and tubing are built-in to withstand vacuum level of up to 500mmHg
- Inner suction cannula has open distal end
- ► 210cm in length with Fr. 24 connector

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-PP-SSSPD5	Simplex Poole Abdominal Drain, 210cm	Ea	1	50	05111





Suction Connecting Tube

Simplex Suction Connecting Tube is used for connecting a suction device to the suction container or the vacuum source

FEATURES AND BENEFITS:

- ► 210cm transparent tubing allows clear visualization of fluid
- Kink-resistant tubing for continuous flow
- ► With integral connectors at both ends for convenient and secure connection
- Universa

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-PP-SCT000	Simplex Suction Connecting Tube	Ea	1	50	MDRE





Suction Connecting Tube

Unomedical Suction Connecting Tube is used as an integration between the suction device and the collecting container.

FEATURES AND BENEFITS:

- ► Made from medical grade PVC
- ► Can withstand high vacuum pressure
- ► Cut-to-fit connector, enables fitting to all types of canisters/suction machines
- ► Connector ensures safe and tight grip
- ► Kink-resistant ensuring unrestricted flow
- ► With 210 cm tubing

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-UMCTFF	Unomedical Suction Connecting Tube F/F	Ea	1	50	04519A
0-PM-UMCTFM	Unomedical Suction Connecting Tube F/M	Ea	1	50	04519A





Universal Funnel Tubing

Medovations Universal Funnel Tubing offers quick connection to tubings or suction units and devices.

FEATURES AND BENEFITS:

- ► Made from clear, kink-resistant medical-grade PVC material that has high-resistance to vacuum collapse
- Available in an easy-dispensing storage
- ► It features a bubble every 3 feet for identification along with quick integral connector
- ► 100ft. tubing

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-PP-MUFT10	Medovations Universal Funnel Tubing	Roll	1	1	MDRE



Pana-Vac®

Tri-Spring Wound Suction Kit with Trocar

Pana-Vac Tri-Spring Wound Suction Kit is a closed wound drainage system, designed to prevent accumulation of fluid which may predispose the wound to infection, which may delay wound healing.

- ► Tri-spring design for easy suctioning
- ► Silicone drain tube with radiopaque line
- ► With universal wound drain y-connector which can be cut to fit to Fr. 7, 10, 15 or 19 drain tube
- ► With anti-reflux valve to prevent fluid backflow
- ► Clear sidewall with graduation for quick and easy measurement observation

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-PP-STS410	Pana-Vac Tri-Spring Wound Suction Kit w/ Trocar, 400cc	Kit	1	20	09635



Pana-Vac®

Closed System Wound Suction Kit

Pana-Vac Closed System Wound Suction Kit is used to remove fluid, blood and exudates from the wound area to prevent infection which may delay wound healing.

FEATURES AND BENEFITS:

- ► Made from silicone material, biocompatible and less traumatic
- ► Injection molded for strength and drain integrity
- ► Barium impregnated for x-ray detection of drain
- ► Clear reservoir and graduations allows easy visualization of measurement
- ► Internal anti-reflux valve, preventing fluid backflow

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available in 100cc, 200cc and 400cc	Set	1	50	01677
Pana-Vac Trocar Fr-15	Ea	1	_	01677



Unimex[®]

Penrose Drain Tubing

Unimex Penrose Drain Tubing is a soft, flexible rubber tube, used as a surgical drain, to prevent the build-up of fluid in a surgical site.

FEATURES AND BENEFITS:

- ► Smooth tubing wall both inside and out
- ► Soft and flexible pure latex rubber material
- ► 18 inches in length
- ► Available in various sizes
- X-ray opaque

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available Diameter Size: 1/4" to 1"	Ea	1	50	1000	08104



simplex

T-Tube Drain

Simplex T-Tube Drain is a soft pliable device with bile stem, designed for post-operative drainage for patients who have undergone gallbladder surgery or surgery of the surrounding tubes draining the gallbladder.

- ► Made of natural latex, silicone coated
- ► 400mm stem length
- ► 200mm crosshead length
- ► Non-toxic, pyrogen-free
- ► Sterilized by ethylene oxide

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-12 to Fr-18	Ea	1	20	200	_



Unimex®

Chest Drainage Thoracic Catheter

Unimex Chest Drainage Thoracic Catheter is designed for easy and atraumatic removal of air or fluids from the pleural cavity to relieve life-threathening chest conditions.

FEATURES AND BENEFITS:

- ► Thermosensitive tube that softens according to the patient's body temperature
- ► X-ray opaque sentinel line, for placement confirmation
- ► Smooth catheter eyes for maximum drainage, while reducing tissue debris build-up
- ► Tapered proximal end, minimizing trauma without tension and compression
- ► Graduated markings every 2cm, ensures accurate placement

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-16 to Fr-36	Ea	25	250	04276



medovations.


Thoracic Catheter

Medovations Thoracic Catheter is designed to maximize patient's comfort with its smooth side-eyes and open distal tip, which allows consistent catheter flow for effective drainage of air or fluids from the pleural space.

FEATURES AND BENEFITS:

- ► Made from clear, kink-resistant, medical grade PVC material
- ► With radiopaque line for placement confirmation
- ► Quick to connect funnel connector provides fast, easy connection to standard drainage units
- ► Available in straight and right angle configurations

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-16 to Fr-36	Pack	1	10	05229





Chest Drainage

Biometrix Chest Drainage System is a compact design chest drainage system, engineered for convenient collection of fluid and air from the pleural space. It features water seal operating system and suction control chamber, which is designed to maintain the required amount of vacuum pressure in the patient's chests.

- ► Clear chambers and tubing for visual drainage observation
- ► With distinct graduation markings for easy measurement reading
- With high negativity float valve, this safety feature prompts the clinicians in the event of high negative pressure
- ► Allows connection of two drainage tubes, preventing the use of Y-connector
- ► Convenient and easy to set-up

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-BM-NC3311	Biometrix Chest Drainage, 1150cc	Ea	1	5	MDRE
2-PN-NC5512	Biometrix Chest Drainage, 2300cc High Cap	Ea	1	5	MDRE





Safeliner

De Royal Safeliner provides simple, yet innovatively safe fluid management for healthcare facilities by allowing safe handling of medical waste fluids, improving infection control practice and protecting hospital investment.

FEATURES AND BENEFITS:

- Innovative adhesive seal on lid keeps the lid securely in place to help prevent leakage and spills once full. Spill-Safe for clinicians
- ► Self-aligning lid allows quick and easy lid assembly.
- ► Improved flow rates allow greater suction power in clinical settings.
- ► Innovative filter technology helps prevent backflow of aspirated fluids into the vacuum system.
- ► Canisters are nested for efficient packaging and storage
- ► Available in 1.5L and 3L
- ► 18" funnel tubing, partner with tandem manifold for automatic flow transfer once one liner is full.
- ► Tandem manifold, to tandem up to four canisters, ideal for areas that requires large volume suctioning.
- ► Canister stand with casters, allows flexibility in location and configuration

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-DRSRL1	Available Size: 1500cc and 3000cc	Ea	1	100	MDRE
1-DR-715046	DeRoyal Tandem Manifold	Ea	1	_	MDRE
1-DR-711142	DeRoyal 18" Funnel Tubing	Ea	1	100	MDRE
1-DR-715045	DeRoyal Stand with Casters	Ea	1	_	MDRE







Vessel Loop

DeRoyal Vessel Loops are used in surgical situation for identification, occlusion and retraction of veins, arteries, ureters and tendons.

FEATURES AND BENEFITS:

- ► Made from high-quality medical-grade silicone
- ► Non-porous to slide smoothly and easily over vessels
- ► Color-coded to assist in vessel detection
- ► Radiopaque X-Ray detectable
- ► Neatly wrapped silicone vessel loops on a plastic card
- ► Minis are 1.3mm wide and 0.9mm thick

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in White, Red, Blue and Yellow	Box	2	10	MDRE





Solidifier Gel

DeRoyal Solidifier is used to quickly and completely solidify blood and body fluids for safe disposal. Solidifies potentially serious infectious waste within the suction canister, quickly and effectively.

- Quickly penetrates throughout the canister to solidify blood and body fluids
- ► Gels in seconds, ten times faster than other solidifiers. No shaking required
- Retains solidity to prevent spills and leaks for safe and easy bagging, transport and disposal
- ► Concentrated formula: one ounce cost effectively solidifies 1000 cc of liquid, twice as much as other products
- ► Chlorine-free for non-toxic incineration

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-DR-711800	DeRoyal Solidifier Gel, 18000cc	Ea	1	12	MDRE



Angiosenz[®]

3 Port Manifold

Angiosenz 3 Port Manifold is a high-pressure manifold, used to manage fluid flow during angiography procedures.It is designed to control fluid flow without disconnecting any existing intravenous line when connecting new lines. It can be connected to saline infusion, contrast medium, and monitoring system during the procedure.

FEATURES AND BENEFITS:

- ► Made from PTFE standard material, tapered and siliconized
- ► Ergonomic design for one hand operation
- No additional fluid lines necessary
 Fewer connections means reducing air bubbles
- ► Clear housing for visual observation
- ► Luer lock fitting male / female connectors
 ► Easy to set up, enhances time management during procedures

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
2-SM-3PMRTO	Angiosenz 3 Port Manifold, Right-On	Ea	40	160	_





2 & 3 Port Manifold

Biometrix 2 & 3 Port Manifold is designed for convenient use during angiography procedures.

FEATURES AND BENEFITS:

- ► Made from PTFE-standard material that is tapered and siliconized
- A flexible, easy-grip wing is provided for ease of placement
- ► Clear housing for visual observation
- Luer lock fitting male / female connectors
- ► Non-pyrogenic

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Av ailable in 2 and 3 Port Manifold	Ea	1	50	_




Gigli Saw Wire

Gigli Saw Wire is a flexible wire saw used by surgeons for bone cutting. Used mainly for amputation, where the bones have to be smoothly cut at the level of amputation.

FEATURES AND BENEFITS:

- ► 9 inches in length
- ► Made from stainless steel material
- ► Premium OR-Grade
- ► Latex-free
- ► Non-sterile

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available Size: 12in to 30in	Ea	1	20	_



THOMAS

Hydrophobic Suction Filter

Thomas Hydrophobic Suction Filter prevents fluid and aerosol contamination from entering and damaging the pump.

- $\begin{tabular}{ll} \begin{tabular}{ll} \beg$
- ► Repels water

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
4-EP-TSMF00	Thomas Hydrophobic Suction Filter 641213	Ea	1	12	_



National Scrub Brush

National Scrub Brush is used for surgical scrubbing prior to any surgical procedure, with resilient long lasting nylon bristles that are resistant to chemicals and provides thorough scrubbing.

FEATURES AND BENEFITS:

- ► Bristles withstand harsh chemicals and autoclavings
- ► Economical, reusable and autoclavable
- ► 112 bristles, 7x16 tufts

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
0-GE-NTLSBH	National Scrub Brush	Ea	1	12	576	MDRE





Surgical Needles

Anchor Surgical Needles are used for threading sutures during surgical operations.

- ► Made from stainless steel material
- ► With comprehensive types and size for different surgical operations
- ► Made in U.S.A

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-AM-1851	Anchor Surgical Needle Abscission 3/8 Rev Cut	Doz	_	_	05148
1-AM-1854	Anchor Surgical Needle Bunnels Str Cut	Doz	_	_	05148
1-AM-1835	Anchor Surgical Needle Cervix 3/8 Rnd	Doz	_	_	05148
1-AM-1830	Anchor Surgical Needle Eye 3/8 Rev Cut	Doz	_	_	05148
1-AM-1831	Anchor Surgical Needle Eye 1/2 Rev Cut	Doz	_	_	05148
1-AM-1833	Anchor Surgical Needle Eye 1/2 Rnd	Doz	_	_	05148
1-AM-1842	Anchor Surgical Needle Ferguson 1/2 Rnd	Doz	_	_	05148
1-AM-185303	Anchor Surgical Needle Fish Hook Rev Cut	Doz	_	_	05148
1-AM-1832	Anchor Surgical Needle Fistula 1/2 Rev Cut	Doz	_	_	05148
1-AM-1861	Anchor Surgical Needle French Eye 1/2 Rnd	Doz	_	_	05148
1-AM-1852	Anchor Surgical Needle Hagedorn Str Cut	Doz	_	_	05148
1-AM-1829	Anchor Surgical Needle Hernia 1/2 Blunt	Doz	_	_	05148
1-AM-1847	Anchor Surgical Needle Hernia Koontz 1/2 Rnd	Doz	_	_	05148
1-AM-1843	Anchor Surgical Needle Intestinal Str Rnd	Doz	_	_	05148
1-AM-1850	Anchor Surgical Needle Intestinal 1/2 Rnd	Doz	_	_	05148
1-AM-1857	Anchor Surgical Needle Jameson Strabismus Cut	Doz	_	_	05148
1-AM-1827	Anchor Surgical Needle Keith Str Trngl	Doz	_	_	05148
1-AM-1838	Anchor Surgical Needle Kidney 1/2 Cvd Blunt	Doz	_	_	05148
1-AM-1839	Anchor Surgical Needle Kidney 1/2 Blunt	Doz	_	_	05148
1-AM-1840	Anchor Surgical Needle Kidney Str Blunt	Doz	_	_	05148
1-AM-1841	Anchor Surgical Needle King 3/8 Rev Cut	Doz	_	_	05148
1-AM-1844	Anchor Surgical Needle Lanes Cleft Palate 1/2 Rev Cut	Doz	_	_	05148
1-AM-1848	Anchor Surgical Needle Lanes Cleft Palate 3/8 Rev Cut	Doz	_	_	05148
1-AM-1846	Anchor Surgical Needle Loopuyt's Rev Cut	Doz	_	_	05148
1-AM-1860	Anchor Surgical Needle Martin Uterine 1/2 Rev Cut	Doz	_	_	05148
1-AM-1823	Anchor Surgical Needle Mayo Intestinal 1/2 Rnd	Doz	_	_	05148
1-AM-1824	Anchor Surgical Needle Mayo 1/2 Rnd	Doz	_	_	05148
1-AM-1828	Anchor Surgical Needle Mayo Intestinal Str Rnd	Doz	_	_	05148
1-AM-1837	Anchor Surgical Needle Intestinal 3/8 Rnd	Doz	_	_	05148
1-AM-1859	Anchor Surgical Needle McGowan 1/2 Rnd	Doz	_	_	05148
1-AM-1825	Anchor Surgical Needle Murphy Intestinal 1/2 Rnd	Doz	_	_	05148
1-AM-1821	Anchor Surgical Needle Regular Surgeon's 1/2 Cvd Rev Cut	Doz	_	_	05148
1-AM-1822	Anchor Surgical Needle Regular Surgeon's 3/8 Rev Cut	Doz	_	_	05148
1-AM-1834	Anchor Surgical Needle Regular Surgeon's 1/2 Rev Cut	Doz	_	_	05148
1-AM-1845	Anchor Surgical Needle Tonsil Davis 3/8 Rnd	Doz	_	_	05148
1-AM-1849	Anchor Surgical Needle Tonsil Davis 1/2 Rnd	Doz	_	_	05148
1-AM-1826	Anchor Surgical Needle Trocar 1/2 Trc	Doz	_	_	05148
1-AM-1836	Anchor Surgical Needle Trocar 3/8 Trc	Doz	_	_	05148






simplex

Nelaton Catheter

Simplex Nelaton Catheter is inserted into the urethra and onto the bladder, designed to allow intermittent bladder catheterization.

FEATURES AND BENEFITS:

- ► Made from non-toxic medical grade PVC
- ► Distal end is coned for atraumatic introduction
- ► Color-coded for easy size identification
- ► 40cm length




DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-8 to Fr-18	Ea	1	50	500	03743



Nelaton Catheter

Unomedical Nelaton Catheter is used for intermittent bladder catheterization.

- ► Color-coded for easy size identification
- ► Made from medical grade PVC material that provides smooth finish for a traumatic introduction
- ► Large drainage lumen provides fast and effective drainage
- ► Transparent silk surface and smoothly rounded tip provides the best possible protection of the sensitive tissue





Urosenz[®]

Red Nelaton Catheter

Urosenz Red Nelaton Catheter are used for intermittent emptying of the bladder for either male or female.

FEATURES AND BENEFITS:

- ► Silicone coated material, preventing use of lubricating gel
- Soft rounded tip, safe to the surrounding tissues
- With 2 lateral eyes ensuring maximum drainage without clogging

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-8 to Fr-18	Box	10	50	06016



Urosenz[®]

Silicone Foley Balloon Catheter

Urosenz Silicone Foley Catheter is made from 100% silicone, biocompatible and completely inert, preventing tissue irritation and encrustation during extended period of indwelling bladder

FEATURES AND BENEFITS:

- ► With radiopaque line and tip which allows visualization during x-ray, for placement confirmation
- ► Clear design of the body and funnel allows to see potential blockage
- ► Large internal lumen for maximum drainage
- Symmetrical balloon which provides strength, symmetry and elasticity
- Needleless inflation port for prevention of needlestick injury



DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-8 to Fr-24	Ea	10	500	05109


Urosenz[®]

Latex Foley Balloon Catheter

Urosenz Latex Foley Balloon Catheter is a flexible tube that pass through the urethra and into the bladder to drain the urine. Ideal for short term to medium term use ranging three to five days.

- ► Made from permanent silicone elastomer that protects against urethral irritation and tissue encrustation
- Exceptionally smooth finish for easier insertion and maximum patient comfort
- Strong catheter shaft with optimal wall thickness that provides strength and flexibility
- Pressure resistant balloon that allows exceptional seal in the patient's bladder

- Needleless inflation valve ensures prevention of needlestick injury
- With additional irrigation port for 3-way Foley Balloon Catheter





simplex

Foley Balloon Catheter (Silicone Coated)

Simplex Foley Balloon Catheter is a flexible tube, designed for convenient effective drainage of urine from the bladder.

FEATURES AND BENEFITS:

- ► Made from silicone elastomer bonding with exceptionally smooth finish for easier insertion & maximum patient comfort.
- ► Distal end is coned for atraumatic introduction
- ► Color-coded for easy size identification
- ► Pressure resistant balloon allows easy inflation or deflation, maximizing patient safety.



DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available Size: Fr-8 to Fr-24	Ea	10	500	05022C

Unimex[®]

Male External Condom Catheter

Unimex Male External Condom Catheter is a flexible latex sheath that slides over the penis and is connected to a urine bag, providing urine passage for male patients. It is a non-invasive catheter, used by men to address urinary incontinence.

FEATURES AND BENEFITS:

- ► Made from soft and lightweight latex material
- ► With hypoallergenic powder, reduces risk for skin irritation
- ► With double sided adhesive foam strap that keeps the condom in place

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Available in Small, Medium and Large	Ea	50	1000	06150



Unimex[®]

Pediatric Urine Collector

Unimex Pediatric Urine Collector is designed for urine collection on infant and pediatric patients.

- ► Comfortable foam padding adhesive that allows secure fit and easy removal from soft skin of infants
- ► Graduated plastic bag for quick urine output measurement
- ► latex-free, durable and easy-to-use
- ► 100ml capacity





Urosenz[®] Urine Leg Bag

Urosenz Urine Leg Bag is a urine collecting bag with two straps, conveniently placed in one leg, designed to provide discreet urine storage for mobile patients.

FFATURES AND BENEFITS:

- Pre-attached 10cm tubing with universal connector that fits all sizes of foley catheters
 Kink-resistant tubing that permits continuous flow of urine
- Flexible garter straps with perforations and buttons that holds the leg bag in place
- Latex-free bag with 800ml capacity, made from transparent PVC material that allows observation of urine color and flow
- One-way valve that prevents backflow of urine to prevent infection
- Easy to operate twist drainage port that provides fast emptying of urine bag

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-UN-ULG800	Urosenz Urine Leg Bag	Ea	1	50	_



ےimplex ـ

Urine Bedside Bag

Simplex Urine Bedside Bag is a simple and reliable urine collection bag, ideal for short term or long

FEATURES AND BENEFITS:

- ► Transparent bag with calibration markings for accurate measurement
- ► With plastic hanger for secure and comfortable placement
- Kink-resistant tubing that allows continuous urine flow

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Available in 2L: 90cm or 120cm Tubing	Ea	10	200	08323



Urosenz[®]

Unigard Urine Drainage Bag with non-return valve

Urosenz Unigard Urine Drainage Bag with non-return valve an innovative type of urine bag that features a needleless sample port, designed for needlefree urine sample taking.

FEATURES AND BENEFITS:

- ► Made from non-toxic, graduated medical-grade plastic
- ► 100% leak-proof and puncture resistant bag with reinforced attachment holes and strap for easy hanging
- ► With adjustable bed sheet clamp
- ► With universal transparent connector
- ► Sample port that ensures needlefree sample taking
- ► Perfect for short-term and long-term urine collection
- ► Slide slip drainage port allows easy, one hand urine emptying
- ► 2500ml bag capacity with 120cm tubing
- ► Equipped wtih non-return valve to prevent retrograde infection

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-SM-UDBA00	Urosenz Unigard Urine Drainage Bag	Ea	1	20	05724



Urosenz°

Urine Meter

Urosenz Urine Meter offers an efficient, high-quality and cost economical closed urine drainage system, designed with convenient to use chamber that provides accurate hourly monitoring of patient's urine output.

- ► Needle-free sample port allows needle-free urine sampling, thus preventing needlestick injury
- ► Built-in non-return valve to effectively prevent retrograde infection
- ► Urine chamber with large scale and easy to read digits provides accurate monitoring of urine output
- ► Easy-to-operate bottom outlet ensures quick and hygienic urine bag emptying

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Urosenz Urine Meter UM1	Ea	1	10	871387





Unometer Safeti Plus

Unomedical Unometer Safeti Plus is a closed urine drainage system, equipped with non-return valve and needleless sample port; used in operating room, intensive care unit, emergency department and other specialty ward for continuous monitoring of urine output in critically ill patients.

- Equipped with non-return valve (NRV) which prevents residual urine from re-entering the foley catheter
- ► Kombikon sample port that is made of rubber material which allows needle-free urine sampling
- Double lumen allows free flow of urine into the collecting bag, which minimizes the risk of urine stasis in the tubing
- ► Easy-to-read scale with large digits for accurate measurement
- ► 500ml chamber capacity and 2000ml bag capacity



REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-PM-510190	Unomedical Unometer Safeti Plus	Ea	1	10	883029

So hurry, book your orders now!



www.inmed.com.ph



Pay by cash, cheque or hussle-free payment options:

Online and over-the-counter payment through banks and credit cards.

INMED CORPORATION

5 Calle Industria, Bagumbayan, Quezon City 1110, Philippines















