

**20<sup>TH</sup>**  
**ANNIVERSARY**  
**CATALOG**  
2023 EDITION



CLEANSPACE®  
POWERED AIR PURIFYING RESPIRATOR (PAPR)



The first ISO 9001 certified medical products trading company.

**INMED**   
CORPORATION

We help you grow your business.



## Who we are

Inmed Corporation is the country's leading and trusted master distributor of quality local and imported medical products and supplies. With the primary focus of providing solid marketing support, Inmed Corporation has established a dependable structure offering significant value to its resellers and dealers in the Philippines and neighboring Asia Pacific Countries.

Committed to help both local and multinational medical partner companies in the Philippines to achieve their target goals. Inmed Corporation continues to work hand-in-hand with them to extend and further expand market reach. The company utilizes its financial strength to purchase high volume supplies and offer its resellers savings at the best attainable cost.

Managed by a team of top professionals in the industry with combined experience of over 60 years in medical products, sales, marketing and distribution, Inmed Corporation boasts of seasoned round-table negotiations, financial packaging and the best business practices benchmarked on international standards providing economical and fundamentally sound solutions to your business needs.

## Our Mission

Our mission is to provide quality local and imported medical products at competitive prices and offer professional customer-oriented service.

## Our Vision

Our vision is to be the Philippines' leading master distributor of quality medical products – being the preferred supplier of all resellers in the industry and recognized as a strong partner in sustainable healthcare in the nation.

## ISO Certified

INMED Corporation is ISO 9000:2015 certified by TUV, compliant with the International Organization for Standardization based in Geneva, Switzerland. It is our commitment to provide our customers with only high quality products, unparalleled service excellence, and reliable timely delivery of all customer orders. Customer satisfaction is our guarantee.



# TABLE OF CONTENTS

Half Mask ver2	8
HALO Cleaning and Storage Plug	6
HALO Half Mask	5
HALO Head Harness	6
HALO Mask ID Badges	7
HALO Particulate Filter	6
HALO Power Unit	5
HALO Unit ID Panels	8
Hood with Visor	7
Particulate Filter P3	9
Particulate Pre-Filter Standard	9
Powered Respirator ver2	8
Protective Coverall and Pre-Filter	6
Steri-Plus Exhalation Valve Filter	7
Steri-Plus Exhalation Valve Filter Case (Reusable)	7



**CleanSpace®**  
RESPIRATORS

## POWERED AIR PURIFYING RESPIRATOR (PAPR)

**NIOSH**  
National Institute for  
Occupational Safety and Health  
CERTIFIED & APPROVED

**CE**

**IP66**

**AS**

**Australian Standard**  
AS/NZS 1716:2012-PAPR-P2  
Issued By BSI  
BMP: 69398



# CleanSpace®



## HALO Power Unit

CleanSpace® HALO is a Powered Air Purifying Respirator (PAPR) system that houses smart technology, designed by biomedical engineers in a revolutionary compact system. A lightweight compact powered respirator with minimal parts (no cables, hoses or helmets), engineered to deliver the highest protection in healthcare settings while providing comfort, easy wearing and easy integration into any settings.

### FEATURES AND BENEFITS:

- ▶ AirSensit™ Technology, a smart automated technology that senses the wearer's breathing to trigger the airflow to ON or OFF.
- ▶ Provides High Protection: HEPA filter provides 99.97% filtration efficiency for particles 0.3µm and larger, including biohazards (APF 50 & 1000).
- ▶ Compact design and lightweight (350g/1lb), compatible with standard hospital PPE. No belts, hoses, battery packs or head tops worn above the shoulders.
- ▶ Comfortable to wear for extended periods, don/doff in under 10 seconds.
- ▶ Easy-to-clean/disinfect and integrates into standard protocols
- ▶ Water tolerant with IP66 rating
- ▶ Longer battery usage: Up to 9 hours of continuous use with rapid charging time in less than two hours
- ▶ Sturdy, shock-resistant casing ensures durability of the unit.
- ▶ Visual/Audible alarms for: particulate filter end of life, battery charge level, low battery warning and manufacturer's minimum design flow.
- ▶ Inclusions: 1 Filter, 1 Head Harness, 1 Fitting Seal Check Cap, Charger and Adaptors
- ▶ Proudly Made in Australia

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-CS3000	Cleanspace Halo Power Unit	Box	1	—	—

# CleanSpace®

## HALO Half Mask

CleanSpace HALO Half Mask is specifically designed to complete the CleanSpace Halo Power Unit, made from soft silicone material, readily fits to the face and does not cause allergic dermatitis.

### FEATURES AND BENEFITS:

- ▶ Made from soft medical grade silicone material, latex free, providing comfortable mold seal.
- ▶ Transparent and facilitates clear voice transmission for easy communication.
- ▶ Easy to connect and remove from the power unit.
- ▶ Compatible with most cleaning and disinfection protocols.
- ▶ Available in 3 sizes to facilitate good seal and comfortable fit.
- ▶ Proudly Made in Australia



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-CS3003	Cleanspace Halo S3 Half Mask, Small	Ea	1	—	—
1-CS-CS3004	Cleanspace Halo S3 Half Mask, Medium	Ea	1	—	—
1-CS-CS3005	Cleanspace Halo S3 Half Mask, Large	Ea	1	—	—

## CleanSpace® HALO Head Harness

CleanSpace HALO Head Harness is used with the CleanSpace Halo PAPR Power Unit and half mask to provide secure, comfortable and convenient fit to the user.

### FEATURES AND BENEFITS:

- ▶ Strap is made from high-density polyethylene (HDPE) material
- ▶ With multiple adjustment points
- ▶ Easy to wash and allows quick drying
- ▶ Durable and designed to provide comfort during long wearing in harsh environments
- ▶ Allows easy and quick replacement

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-CS3008	Cleanspace Halo Head Harness	Ea	1	—	—



## CleanSpace® HALO Particulate Filter

CleanSpace HALO Particulate Filter P3 provides suitable protection against airborne particulates (dust, mists and fumes) and other biohazards.

### FEATURES AND BENEFITS:

- ▶ Used with the revolutionary CleanSpace PAPR
- ▶ HEPA/P3 protection provides 99.97% filtration efficiency for 0.3micron and larger.
- ▶ Suitable for protection against airborne particulate (dust, mists and fumes).
- ▶ Materials: Fiberglass particulate media, plastic casing and silicone seal
- ▶ Allows easy and quick fitting and removal from the power unit.
- ▶ These filters are not waterproof and should be replaced if they come in contact with water.
- ▶ Pack of 3
- ▶ Proudly Made in Australia

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-CS3002	Cleanspace Halo Particulate Filter P3 Standard	Pack	3	—	—



## CleanSpace® HALO Cleaning and Storage Plug

CleanSpace HALO Cleaning and Storage Plug is designed to be used with the CleanSpace HALO Respirator Unit to seal all the openings during cleaning and storage.

### FEATURES AND BENEFITS:

- ▶ Used with the revolutionary CleanSpace HALO PAPR
- ▶ Made from polycarbonate and silicone material
- ▶ Weight: 85g
- ▶ Dimensions: 180mm x 110mm x 800mm
- ▶ Proudly Made in Australia

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-CS3011	CleanSpace Halo Cleaning and Storage Plug	Ea	1	—	—



## CleanSpace® Protective Coverall and Pre-Filter

CleanSpace Protective Coverall and Pre-Filter is designed to protect the CleanSpace Respirator and the Particulate P3 Filter against airborne particulates like dust, mists and fumes to enhance the life of the Particulate P3 Filter.

### FEATURES AND BENEFITS:

- ▶ Compatible with any CleanSpace Respirator
- ▶ Made from spun polymer fibers
- ▶ Easy to insert and to remove to any CleanSpace Respirator
- ▶ Sold in pack of 20's
- ▶ Approved Standards: AS/NZS1716: 2012, EN 12942
- ▶ Proudly Made in Australia

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-PF0058	CleanSpace Protective Coverall and Pre-Filter	Pack	20	—	—



## CleanSpace® Reusable Steri-Plus Exhalation Valve Filter Case

CleanSpace Steri-Plus Exhalation Valve Filter Case houses the Steri-Plus Exhalation Valve Filter; together they are designed to filter the exhaled air of the wearer to protect the environment from airborne contaminants.

### FEATURES AND BENEFITS:

- ▶ Used with the revolutionary CleanSpace HALO
- ▶ Materials: Inner Case – ABS, Outer Case – Polycarbonate.
- ▶ Reusable and easy to clean
- ▶ Proudly Made in Australia

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-CS3038	CleanSpace Steri-Plus Exhalation Valve Filter Case	Ea	1	—	—



## CleanSpace® Steri-Plus Exhalation Valve Filter

CleanSpace Steri-Plus Exhalation Valve Filter is suitable for use with the CleanSpace Halo Power System; it is designed to filter the exhaled air of the wearer to protect the surrounding people and the environment from airborne contaminants.

### FEATURES AND BENEFITS:

- ▶ Used with the revolutionary CleanSpace HALO
- ▶ Made from non-woven electrostatic media, polypropylene material
- ▶ Allows quick fitting and removal from the exhalation valve filter case
- ▶ Used alongside with the CleanSpace Particulate Filter
- ▶ Proudly Made in Australia

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-CS3039	CleanSpace Steri-Plus Exhalation Valve Filter	Pack	20	—	—



## CleanSpace® Hood with Visor

CleanSpace Hood with Visor is suitable for use with the CleanSpace Half Mask, it is designed to provide protection against splashes and airborne contaminants.

### FEATURES AND BENEFITS:

- ▶ Used with the revolutionary CleanSpace Respirator
- ▶ Materials: Fabric – microporous polyethylene film laminated to spunbond polypropylene, Visor – PVC.
- ▶ Wide panoramic visor.
- ▶ Proudly Made in Australia

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-PF0097	CleanSpace Hood with Visor, for Half Mask	Ea	1	—	—



## CleanSpace® HALO Mask ID Badges

The CleanSpace Mask ID Badges are suitable for use with the CleanSpace HALO Half Mask. They allow the user to individualize their Half Masks through the addition of a colored badge.

### FEATURES AND BENEFITS:

- ▶ Used with the revolutionary CleanSpace HALO
- ▶ Material: ABS.
- ▶ Allows easy fitting and removal from the CleanSpace HALO Respirator.
- ▶ Three different colors available.
- ▶ Proudly Made in Australia

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-OR3015	CleanSpace Halo Mask ID Badge, Orange (Order in 10s)	Pair	2	—	—
1-CS-YE3015	CleanSpace Halo Mask ID Badge, Yellow (Order in 10s)	Pair	2	—	—
1-CS-BL3015	CleanSpace Halo Mask ID Badge, Blue (Order in 10s)	Pair	2	—	—



## CleanSpace®

### HALO Unit ID Panels

The CleanSpace HALO Unit ID Panels are suitable for use with the CleanSpace HALO. They allow the user to individualize their units by replacing the standard HALO panels with a range of colored ones.

#### FEATURES AND BENEFITS:

- Used with the revolutionary CleanSpace HALO
- Material: PC/ABS.
- Allows easy fit and removal from the CleanSpace Halo
- Proudly Made in Australia



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-BL3013	CleanSpace Halo Unit ID Panels, Blue	Pair	2	—	—
1-CS-GR3013	CleanSpace Halo Unit ID Panels, Green	Pair	2	—	—
1-CS-PI3013	CleanSpace Halo Unit ID Panels, Pink	Pair	2	—	—
1-CS-PU3013	CleanSpace Halo Unit ID Panels, Purple	Pair	2	—	—
1-CS-DP3013	CleanSpace Halo Unit ID Panels, Dark Pink	Pair	2	—	—
1-CS-OR3013	CleanSpace Halo Unit ID Panels, Orange	Pair	2	—	—
1-CS-SI3013	CleanSpace Halo Unit ID Panels, Silver	Pair	2	—	—
1-CS-YE3013	CleanSpace Halo Unit ID Panels, Yellow	Pair	2	—	—

## CleanSpace2

### Powered Respirator

The CleanSpace2 Powered Respirator is a battery-powered, positive pressure air-purifying/filtering respirator. This respirator is breathe responsive and is designed to provide comfortable respiratory protection against gases, vapors, bacteria, viruses and other airborne biohazards.

#### FEATURES AND BENEFITS:

- Lightweight, comfortable and easy to wear over long period.
  - Portable: no hoses, belts or waist mounted battery packs.
  - Uses High Efficiency Particulate Filters (Fiberglass Pleated Filter); suitable for use against solid and liquid particulates including dusts, mists, fumes, fibers, bacteria, viruses and other biohazards
  - AirSensit™ System for mask pressure control and breathe responsive airflow.
  - One-button system, standby and auto-start/stop modes for quick don/doffing.
  - Easy-to-use head harness support
  - Long operating time: up to 8 hours and quick battery recharging: <2hours
  - Visual/Audible alarms for: particulate filter end of life, battery charge level, low battery warning, manufacturers minimum design flow (120L/min).
- Proudly Made in Australia



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-LI-CS20034	CleanSpace2 Powered Respirator	Box	1	—	—

## CleanSpace2

### Half Mask

CleanSpace2 Half Mask is specifically designed to complete the CleanSpace2 PAPR power system, made from soft silicone material that readily fits to the face and does not cause allergic dermatitis to the wearer.

#### FEATURES AND BENEFITS:

- Compatible with CleanSpace2 Powered Respirator
- Available in 3 sizes to facilitate good seal and comfortable fit
- Materials: Soft silicone cushion, nylon clips and polycarbonate (exhalation valve)
- Transparent and facilitates communication.
- Lightweight. Weight: average 125g
- Proudly Made in Australia



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-LI-C20033	CleanSpace2 Half Mask, Small	Ea	1	—	—
1-LI-C21010	CleanSpace2 Half Mask, Medium	Ea	1	—	—
1-LI-C20027	CleanSpace2 Half Mask, Large	Ea	1	—	—



## CleanSpace2

### Particulate Filter P3

CleanSpace Particulate Filter P3 is suitable for protection against airborne particulates such as dust, mist and fumes.

#### FEATURES AND BENEFITS:

- Used with the revolutionary CleanSpace2 PAPR
- Materials: Fibreglass particulate media, plastic casing, silicone seal
- Weight: average 21g. Dimensions: 120mm x 45mm x 45mm
- Quickly fitted and easy to put and remove from power unit
- Suitable for protection against airborne particulate (dust, mists and fumes).
- Easy and quickly fitted and removed from the power unit.
- Pack of 3's
- Proudly Made in Australia



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-CS-CS3038	CleanSpace Steri-Plus Exhalation Valve Filter Case	Ea	1	—	—

## CleanSpace2

### Particulate Pre-Filter Standard

CleanSpace Pre-Filter is suitable for use with the CleanSpace Standard Particulate P3 Filter for additional protection against airborne particulate (dust, mist and fumes). This Pre-filter is designed to remove the course particles and enhance the life of the Particulate P3 filter.

#### FEATURES AND BENEFITS:

- Used with the revolutionary CleanSpace2 Power Unit
- Easy to put and removed from the power unit
- Material: Spun polymer fiber
- Sold in pack of 10 pre-filters, intended to be disposed every after use.
- Proudly Made in Australia



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-LI-C20036	CleanSpace Particulate Pre-Filter Standard	Pack	10	—	—

So hurry, book your orders now!



[www.inmed.com.ph](http://www.inmed.com.ph)



**Pay by cash, cheque or hassle-free payment options:**

Online and over-the-counter payment through banks and credit cards.

## INMED CORPORATION

5 Calle Industria, Bagumbayan, Quezon City 1110, Philippines

☎ 02.85711888 📞 📞 **Globe** 0917.8146633

Follow us

